

RANDOX

RANDOX INTERNATIONAL
QUALITY ASSESSMENT SCHEME



RIQAS

RANDOX

QUALITY CONTROL

RIQAS

THE LARGEST INTERNATIONAL EQA SCHEME
WITH OVER 50,000 LAB PARTICIPANTS



| | |
|----|--|
| 01 | BENEFITS |
| 02 | EQA |
| 03 | RIQAS REPORTS |
| 04 | WEB-BASED DATA TRANSFER |
| 05 | PARTICIPATION IN RIQAS |
| 06 | STANDARD REPORT |
| 07 | TEXT SECTION |
| 08 | HISTOGRAM |
| 09 | LEVEY-JENNINGS CHART |
| 10 | TARGET SCORE CHART |
| 11 | %DEVIATION CHARTS |
| 12 | MULTI-METHOD STAT SECTION |
| 13 | SUMMARY PAGE |
| 14 | MULTI-INSTRUMENT REPORT |
| 15 | URINE TOXICOLOGY REPORT |
| 16 | URINE TOXICOLOGY REPORT SCREENING SECTION |
| 17 | URINE TOXICOLOGY REPORT QUANTITATIVE SECTION |
| 18 | URINALYSIS REPORT |
| 19 | SEROLOGY: SCREENING (QUALITATIVE) REPORT |
| 20 | SEROLOGY: SCREENING (QUANTITATIVE) REPORT |
| 21 | QUANTITATIVE END-OF-CYCLE REPORT |
| 22 | CHART SECTION (END-OF-CYCLE REPORT) |
| 23 | TEXT SECTION (END-OF-CYCLE REPORT) |
| 25 | CURRENT & PREVIOUS CYCLE ABSOLUTE SDIs (END-OF-CYCLE REPORT) |
| 26 | CERTIFICATE OF PERFORMANCE (END-OF-CYCLE REPORT) |
| 27 | MONITORING EQA PERFORMANCE |
| 30 | RIQAS PROGRAMMES |
| 35 | PARAMETER INDEX |
| 49 | RELATED PRODUCTS |
| 51 | RANDOX - A GLOBAL DIAGNOSTIC SOLUTIONS PROVIDER |
| 52 | CONTACT US |

BENEFITS

Delivering a comprehensive yet cost effective EQA solution, RIQAS will help meet regulatory requirements and increase confidence in test system accuracy.



Large Database of Users

- A high level of participation means peer group numbers are maximised whilst ensuring availability of data for a wide range of instruments and methods.



User-friendly Reports

- Simple, one page per parameter format, enables at-a-glance performance assessment, saving valuable laboratory time.
- Complimentary multi-instrument and interlaboratory reports allow comparative performance assessment of all laboratory systems and multiple connected laboratories.
- End-of-Cycle reports, summarising performance compared to the previous cycle, allow you to identify improvements in quality over time.



Cost Effective

- Our extensive range of multi-analyte programmes will reduce the number of individual programmes required to cover your test menu, saving both time and money.
- Reduced parameter options for selected programmes offer greater flexibility, ensuring suitability for laboratories of all sizes and budgets.
- Register up to five instruments per programme (volume permitting) at no extra cost for comparative performance assessment.



Frequency

- Frequent reporting allows early identification of system errors and implementation of any necessary corrective actions with minimum disruption to the lab.
- With a turnaround of less than 72 hours for most reports, corrective action can be implemented earlier, potentially reducing costly errors with patient results.



High Quality Samples

- Samples spanning clinically relevant levels allow identification of concentration related biases, helping to ensure accurate instrument performance.
- Human samples free from interfering preservatives increase confidence that EQA performance mirrors the performance of patient samples.
- Reference method values are provided in the Clinical Chemistry programme for selected parameters and lots, while for the Immunosuppressant programme they are provided for all parameters and lots.




Highly Accredited

- Programmes accepted by National and International accreditation bodies worldwide.
- Participant certificates provide evidence of participation in a reputable EQA scheme.

RIQAS is the largest international EQA scheme in the world. It is used by more than 50,000 laboratory participants in 139 countries. 33 programmes are currently available.

RIQAS Programmes

- Ammonia/Ethanol
- Anti-TSH Receptor
- Blood Gas
- BNP
- Cardiac
- Cardiac Plus ^{*coming in 2021}
- Cerebrospinal Fluid (CSF)
- Clinical Chemistry
- Coagulation
- CO-Oximetry
- CYFRA 21-I
- ESR
- Glycated Haemoglobin (HbA1c)
- Haematology
- Human Urine
- Immunoassay
- Immunoassay Speciality I
- Immunoassay Speciality 2
- Immunosuppressant Drugs
- Lipid
- Maternal Screening
- Serology Epstein Barr Virus (EBV)
- Serology (HIV/Hepatitis)
- Serology (Syphilis)
- Serology (ToRCH)
- Specific Proteins
- Sweat Testing
- Therapeutic Drugs
- Trace Elements in Blood
- Trace Elements in Serum
- Trace Elements in Urine
- Urinalysis
- Urine Toxicology



Accreditation

- RIQAS provides certificates as proof of EQA participation and performance for laboratory accreditation purposes.
- RIQAS is a UKAS accredited Proficiency Testing Provider, No. 0010, and is accredited to ISO/IEC 17043:2010, 'Conformity Assessment- General Requirements for Proficiency Testing'.
- Accreditation to ISO/IEC 17043:2010 highlights the superior quality and excellence of RIQAS.

UK Performance Surveillance

- Recognised by the Joint Working Group on Quality Assurance (JWG QA).
- Recognised by various National Quality Assurance Advisory Panels (NQAAP).

Independent Advisory Panel

RIQAS participants have access to an independent advisory panel consisting of scientific and clinical experts. This ensures professional and ethical conduct of the scheme and participant confidentiality.


 RIQAS support staff are on hand to offer
 advice and troubleshoot technical queries.

RIQAS REPORTS

RIQAS reports are presented in a user-friendly, one page per parameter format. This allows easy interpretation of your analytical performance.

RIQAS Reports

- Statistical breakdown by all methods, your method and, where applicable, your instrument, including running means for the last 10 samples.
- Compare your instrument group, method group and all methods using the histogram.
- Identify trends, biases and precision problems using the visual charts.
- The Target Score chart grades your performance in a moving window over the last 20 samples, including the previous cycle.
- At-a-glance summary page for all parameters in the programme.
- Compare your result with statistically robust consensus means.
- Identify acceptable and poor performance using fit-for-purpose performance indicators:
 - SDI
 - %Deviation
 - Target Score



Summary CSV Files

It is possible to receive an additional summary of your report statistics, acceptable limits and performance indicators as a .csv file for every sample.

Multi-Instrument Reports

Laboratories can register up to five instruments at no extra cost. Individual reports for each instrument plus a unique multi-instrument report are provided. The multi-instrument report plots the performance of each individual instrument on a single, colour coded Levey-Jennings chart, ensuring instant identification of any differences in instrument performance. Additional sample packs may be ordered as required if volume supplied is insufficient for the registered instruments.





Laboratory Group Reports

The Group Reporting facility enables laboratory groups or chains to monitor the performance of satellite sites. Each affiliated laboratory will receive their individual reports with the group supervisor also receiving a summary report comparing each laboratory in the network.

WEB-BASED DATA TRANSFER

RIQAS.Net offers easy, direct access for the submission of results and retrieval of reports direct from the RIQAS host server.

- Available in multiple languages.
- Confidentiality and security is maintained through the use of password protected access.
- Submit current, corrected and future results (normal policies apply), directly into the RIQAS database. Receipt of results is confirmed by e-mail.
- Multi-lingual registration identifier provides simple identification of multiple registrations.
- Additions and changes to assay details can be made quickly and easily online.
- Requests for new method, instrument and reagent codes can be made online.
- Reports are emailed in PDF format as soon as they are prepared.
- Reports for the previous two cycles can be downloaded from the website.
- View, print, store or distribute reports as you wish.
- Update your laboratory's certificate of participation details in multiple languages.
- All that is required is web access, Adobe Reader (for viewing reports) and a valid password to access the system.
- No additional software required.




PARTICIPATION IN RIQAS

Participation in RIQAS follows these simple steps:



Performance data is presented in a one page format with up to seven sub-reports.



- 1 **Text Section Chart:** Statistics for all methods, your method and instrument group (programme specific).
- 2 **Histogram Chart:** Method and instrument comparison.
- 3 **Multi-Method Stat Section Chart:** Enables assessment of the performance of each method.
- 4 **Levey-Jennings Chart:** Details features of your laboratory's performance.
- 5 **Target Score Chart:** This unique chart provides a numerical index of performance, allowing at-a-glance assessment.
- 6 **%Deviation by Sample Chart:** Helps to identify trends and shifts in performance.
- 7 **%Deviation by Concentration Chart:** Rapid assessment of concentration related biases.

TEXT SECTION

The text section summarises the statistical information for each parameter.

| Glucose, mg/dl | | 2 | 3 | 4 | 5 | 6 | 7 |
|--|--|---------|---------|-----|----------------|------|------|
| | | N | Mean | CV% | U _m | SDPA | Exc. |
| <input type="checkbox"/> All Methods | | 6598 | 108.157 | 3.9 | 0.06 | 5.33 | 621 |
| <input checked="" type="checkbox"/> Hexokinase | | 3524 | 108.234 | 2.9 | 0.07 | 5.33 | 285 |
| <input checked="" type="checkbox"/> Roche Cobas c501/502 e601/602 e801 | | 677 | 107.618 | 2.2 | 0.12 | 5.30 | 48 |
| ▲ Your Result | | 110.300 | SDI | | 0.51 | | 9 |
| | | | RMSDI | | 0.40 | | 10 |
| ■ Mean for Comparison | | 107.618 | TS | | 101 | | 11 |
| | | | RMTS | | 107 | | 10 |
| 15 Reference Value | | 105.179 | %DEV | | 2.5 | | 12 |
| | | | RM%DEV | | 2.0 | | 10 |
| Acceptable limits derived from Biological Variation | | | | | 6.96% | | 13 |
| Acceptable limits of performance for RIQAS | | | | | 8.10% | | 14 |
| Performance statement appears here if performance indicators exceed limits | | | | | | | |

RIQAS performance indicators include SDI, Target Score and %Deviation.

Acceptable performance criteria:

SDI < 2
Target score ≥ 50
%Deviation < defined acceptable limits

- 1 Report is presented in your chosen unit.
- 2 Number of returned results used to generate Mean for Comparison.
- 3 Average value of all laboratories' results.
- 4 Coefficient of Variation.
- 5 Uncertainty associated with the Mean for Comparison.

$$U_m = \frac{1.25 \times SD}{\sqrt{n}}$$
- 6 SDPA = Standard Deviation for Performance Assessment, calculated from the Target Deviation for Performance Assessment (TDPA) and the Mean for Comparison.

$$SDPA = \frac{TDPA \times \text{Mean for Comparison}}{t\text{-value} \times 100}$$

t-value = factor which represents the % of poor performers reflected in the TDPA (t-value ~ 1.645 when ~10% laboratories achieve poor performance), SDPA is combined with U_m, where appropriate.

If U_m > (0.3 × SDPA) then SDPA_{adjusted} = √(U_m² + SDPA²) and the reported value is suffixed with "a"

If U_m is less than (0.3 × SDPA) then SDPA_{adjusted} = SDPA
- 7 After statistical reduction, some results are excluded from the mean for comparison.
- 8 Ideally this will be your instrument group mean. If N<5 for instrument group, your method group mean is selected as mean for Comparison.
- 9 Standard Deviation Index = $\frac{\text{Your Result} - \text{Mean for Comparison}}{SDPA_{adjusted}}$
- 10 Running Mean average of the last 10 performance indicators is used to monitor performance over time and concentration range.
- 11 Target Score - The closer a value is to 120, the better the performance.

$$TS = \log_{10} \left(\frac{3.16 \times TDPA}{|\%Dev|} \right) \times 100$$
- 12 %Deviation from the Mean for Comparison -


$$\%Dev = \frac{\text{Your Result} - \text{Mean for Comparison}}{\text{Mean for Comparison}} \times 100$$

The closer the value is to zero, the better the performance.
- 13 Biological Variation stated for information purposes only.
- 14 Performance limit set for this parameter.
- 15 Reference values quoted for information purposes, where applicable.

HISTOGRAM

The Bar Graph is intended as a quick visualisation of how your lab's result compares to the method mean, instrument mean and all method mean.

- All methods
- Your method group
- Your instrument group
(programme specific)




- 1** Total of 1126 laboratories reported values between 106.55 and 108.15.

4 Your result is indicated by the black triangle.
- 2** 200 laboratories reported values between 101.77 and 103.36 in your method group.

5 41 laboratories reported values between 111.35 and 112.94 in your instrument group.
- 3** RIQAS reports show your unit of measurement.

LEVEY-JENNINGS CHART

SDIs reflect laboratory performance in relation to fit-for-purpose SDPAs and are useful to monitor performance over time. Acceptable performance is $SDI < 2$.



- 1 The Mean for Comparison for each sample is indicated at the top of the chart. This allows easy assessment of concentration related bias:
 I: Instrument mean
 M: Method mean
 A: All method mean

2 This line indicates a change in registration details for this parameter:

3 Your SDI (Standard Deviation Index).

4 Sample number:


5 N = No result returned from your laboratory.

6 C = Corrected results will be accepted for non-analytical errors. Corrected results will be accepted up to 4 weeks after the final submission deadline, on application, with evidence of analysis. Late results are only accepted if there has been a Randox error.

R = Incorrect results can be removed retrospectively on request.

TARGET SCORE CHART

The Target Score (TS) allows you to assess your performance at a glance. The TS relates the %Deviation of your result from the Mean to a Target Deviation for Performance Assessment (TDPA). TDPAs are set to encourage participants to achieve and maintain acceptable performance. TDPAs are fit-for-purpose performance criteria which are set taking guidance from ISO/IEC 17043, ISO 13528 and IUPAC. Target Deviations for Performance Assessment are also used to calculate the Standard Deviation for Performance Assessment (SDPA).



1 This is the upper deviation limit of performance for this parameter. TDPAs are reviewed regularly and deemed fit for purpose by the RIQAS Advisory Panel.


2 High scores ≥ 50 in the lighter shaded area represent acceptable, good or excellent performance.

3 Heavy shading for values 10 to 50 signifies poor performance.

%DEVIATION CHARTS

The %Deviation by sample chart helps to identify trends and shifts in performance.

$$\%Deviation = \frac{\text{Your Result} - \text{Consensus Mean}}{\text{Consensus Mean}} \times 100\%$$




1 %Deviation from Mean for Comparison.

2 Plot of Running Mean %Deviations (average of the last 10 %Deviations for the sample indicated).

3 Acceptable limits of performance. These are defaulted to RIQAS TDPAs but can be set to e.g. biological variation or regulatory requirement on request.

The %Deviation by concentration chart enables rapid assessment of concentration related biases. Biases at low or high concentrations can be easily determined.



1 Current sample indicated by square.

2 %Deviation at specific concentration.

MULTI-METHOD STAT SECTION

This section provides an easy way of assessing the performance of other methods used to analyse the parameter in question.

| Method | N | Mean | CV% | U_m |
|---------------------------------|----------|-------------|------------|----------------------|
| Hexokinase | 3524 | 108.234 | 2.9 | 0.07 |
| Glucose oxidase | 2673 | 108.474 | 5.5 | 0.14 |
| Ortho Vitros MicroSlide Systems | 251 | 105.036 | 2.3 | 0.19 |
| Glucose dehydrogenase | 63 | 107.832 | 3.5 | 0.59 |
| GOD/02-Beckman method | 37 | 106.425 | 2.5 | 0.55 |
| Other Dry Chemistry | 27 | 108.360 | 6.3 | 1.64 |
| Agappe - GOD-PAP | 21 | 109.727 | 3.7 | 1.11 |
| Oxygen electrode | 17 | 106.666 | 1.3 | 0.43 |
| Vitros, DT60/DT60 II | 3 | 105.595 | 3.9 | 2.99 |

SUMMARY PAGE

Located at the back of the RIQAS Report, the Summary Page collates the key information, allowing participants to review the performance of all parameters at-a-glance.

| Analyte | Mean for Comparison | Your Result | SDI | RMSDI | %DEV | RM%DEV | TS | RMTS | Performance |
|----------------------|---------------------|-------------|--------------|---------|--------------|---------|-----------|------|-------------|
| Albumin | 2.120 | 2.230 | 1.00 | 0.37 | 5.2 | 2.0 | 72 | 107 | |
| Alkaline Phosphatase | 17.705 | 19.000 | 0.61 | -0.27 | 7.3 | -2.9 | 93 | 105 | |
| ALT (GPT) | 12.387 | 12.000 | -0.33 | -0.47 | -3.1 | -3.8 | 119 | 103 | |
| Amylase, Total | 20.454 | 22.000 | 0.72 | -0.29 | 7.6 | -2.5 | 86 | 103 | |
| AST (GOT) | 11.976 | 11.000 | -0.86 | -0.03 | -8.2 | -0.4 | 78 | 100 | |
| Bicarbonate | 8.203 | 6.900 | -1.48 | 0.15 | -15.9 | 1.5 | 54 | 98 | |
| Bilirubin, Direct | 0.251 | 0.380 | <u>2.57</u> | 2.64 | <u>51.3</u> | 47.2 | <u>31</u> | 29 | ▲ |
| Bilirubin, Total | 0.701 | 0.640 | -0.91 | -0.29 | -8.8 | -2.9 | 76 | 101 | |
| Calcium | 6.074 | 6.020 | -0.19 | -0.40 | -0.9 | -1.8 | 120 | 92 | |
| Chloride | 76.353 | 77.000 | 0.30 | -0.28 | 0.8 | -0.8 | 120 | 98 | |
| Cholesterol | 112.696 | 110.000 | -0.55 | 0.05 | <u>2.4</u> | 0.2 | 97 | 115 | |
| CK, Total | 111.659 | 111.000 | -0.08 | 0.35 | -0.6 | 2.5 | 120 | 107 | |
| Creatinine | 0.607 | 0.620 | 0.27 | 0.06 | 2.1 | 0.5 | 120 | 117 | |
| Glucose | 36.429 | 36.000 | -0.26 | -0.84 | -1.2 | -3.7 | 120 | 82 | |
| HDL-Cholesterol | 98.836 | 102.000 | 0.21 | -0.04 | 3.2 | -0.4 | 120 | 113 | |
| Iron | 97.374 | 99.000 | 0.28 | 0.01 | 1.7 | 0.1 | 120 | 114 | |
| Lactate | | No Result | | Too Few | | Too Few | N/A | N/A | |
| LD (LDH) | 85.894 | 87.000 | 0.11 | -0.70 | 1.3 | -6.3 | 120 | 89 | |
| Magnesium | 1.313 | 1.390 | 0.79 | -0.07 | 5.8 | -0.5 | 82 | 107 | |
| Phosphate, Inorganic | 1.451 | 1.540 | 1.02 | 0.02 | 6.1 | 0.1 | 71 | 112 | |
| Potassium | 1.770 | 1.840 | 1.10 | -0.25 | 3.9 | -0.7 | 67 | 99 | |
| Protein, Total | 3.850 | 3.830 | -0.11 | 0.07 | -0.5 | 0.3 | 120 | 114 | |
| Sodium | 112.537 | 114.000 | 0.58 | -0.01 | 1.3 | -0.0 | 95 | 104 | |
| TIBC | 133.143 | 133.000 | -0.01 | -0.01 | -0.1 | -0.1 | 120 | 117 | |
| Trig Total | 23.626 | 24.000 | 0.18 | -0.09 | 1.6 | -0.6 | 120 | 114 | |
| Urea | 5.872 | 5.000 | <u>-2.02</u> | -0.57 | <u>-14.9</u> | -4.0 | <u>41</u> | 95 | ▲ |
| Uric Acid (Urate) | 3.135 | 3.100 | -0.20 | -0.44 | -1.1 | -2.4 | 120 | 107 | |

ORMSDI -0.05

ORM%DEV 0.8


ORMTS 102

- 1 RMSDI - is the Running Mean of the 10 previous SDIs (if fewer than 10 results on file, "Too Few" is printed).
- 2 RM %DEV - Average of the last 10 %DEV for this parameter.
- 3 RMTS - Average of the last 10 Target Scores for this parameter.
- 4 Red triangle appears when all performance indicators (SDI, %DEV and TS) exceed acceptable performance, i.e: when
SDI > 2
TS < 50
%DEV > acceptable limits set

- 5 Overall RMSDI = average RMSDI for this sample distribution.
- 6 Overall RM%DEV = average RM%DEV for this sample distribution.
- 7 Overall RMTS = average RMTS for this sample distribution.

MULTI-INSTRUMENT REPORT

Register up to five instruments per programme at no extra cost. In addition to a standard report for each instrument, a multi-instrument report is also provided allowing comparative performance assessment.



- 1 Allowable %deviation for the parameter in question, based on the RIQAS TDPA.
- 2 Sample number.
- 3 Lab name.
- 4 Unique instrument ID.

- 5 Poor performance.
- 6 %Deviation for each individual sample.
- 7 RM %Dev - Average of the last 10 %Dev for this parameter.
- 8 %Deviation chart comparing the performance of each instrument.

URINE TOXICOLOGY REPORT

Laboratory performance is presented in both quantitative and qualitative screening formats, allowing for easy interpretation at-a-glance.


Screening Section

Quantitative Section

Amphetamines Group, ng/ml

Your Result Negative

Based on weighed-in value of 375
and your chosen cut-off value of 500
the correct response was Negative




All Methods
 KIMS

| N | Mean | CV% | U _m | SDPA | Exc. |
|----|---------|------|----------------|--------|------|
| 45 | 386.683 | 19.1 | 13.78 | 73.94 | 4 |
| 17 | 357.294 | 18.6 | 20.14 | 69.41a | 1 |

| | | | |
|-----------------------|---------|-------|---------|
| ▲ Your Result | 352.000 | SDI | -0.08 |
| | | RMSDI | Too Few |
| ■ Mean for Comparison | 357.294 | | |

| | | |
|---------------|------|-------|
| d-Amphetamine | 375 | ng/ml |
| Ethanol | 45 | mg/dl |
| LSD | 1.25 | ng/ml |
| EDDP | 225 | ng/ml |
| Buprenorphine | 7.5 | ng/ml |

| | Cut-off | TN | TP | FN | FP | RC | NT | Total |
|--|------------|----------|----------|----------|----------|----------|----------|----------|
| Your Result | 500 | 1 | 0 | 0 | 0 | 0 | 0 | 1 |
| KIMS | 300 | 0 | 1 | 0 | 0 | 0 | 0 | 1 |
| | 500 | 12 | 0 | 0 | 0 | 0 | 0 | 12 |
| | 1000 | 9 | 0 | 0 | 0 | 0 | 0 | 9 |
| | All | 21 | 1 | 0 | 0 | 0 | 0 | 22 |
| All Methods | 150 | 0 | 1 | 0 | 0 | 0 | 0 | 1 |
| | 300 | 0 | 7 | 0 | 0 | 0 | 0 | 7 |
| | 500 | 31 | 0 | 0 | 3 | 0 | 0 | 34 |
| | 1000 | 62 | 0 | 0 | 6 | 0 | 0 | 68 |
| | All | 93 | 8 | 0 | 9 | 0 | 0 | 110 |
| Competitive Antibody Binding | 500 | 0 | 0 | 0 | 1 | 0 | 0 | 1 |
| CEDIA | 500 | 4 | 0 | 0 | 0 | 0 | 0 | 4 |
| Chemiluminescence | 500 | 1 | 0 | 0 | 0 | 0 | 0 | 1 |
| DRI-EIA | 500 | 4 | 0 | 0 | 1 | 0 | 0 | 5 |
| EMIT | 500 | 8 | 0 | 0 | 0 | 0 | 0 | 8 |
| FPIA | 500 | 1 | 0 | 0 | 0 | 0 | 0 | 1 |
| Point of Care | 500 | 1 | 0 | 0 | 0 | 0 | 0 | 1 |
| Randox Biochip Array Technology | 500 | 0 | 0 | 0 | 1 | 0 | 0 | 1 |





| Method | N | Mean | CV% | U _m |
|------------------------------|----|---------|------|----------------|
| KIMS | 17 | 357.294 | 18.6 | 20.14 |
| EMIT | 9 | 344.440 | 12.4 | 17.74 |
| Competitive Antibody Binding | 7 | 460.857 | 8.9 | 19.30 |
| DRI-EIA | 7 | 433.767 | 17.2 | 35.31 |
| CEDIA | 2 | 410.000 | 11.7 | 42.25 |

Performance History

| Your Data (Last 10 Samples) | Your Method (This Sample) | Your Method (Last 10 Samples) | All Methods (This Sample) | All Methods (Last 10 Samples) |
|-----------------------------|---------------------------|-------------------------------|---------------------------|-------------------------------|
| 0 % False Negatives | 0 % False Negatives | 1 % False Negatives | 0 % False Negatives | 8 % False Negatives |
| 10 % False Positives | 0 % False Positives | 11 % False Positives | 8 % False Positives | 7 % False Positives |
| 90 % Correct Responses | 100 % Correct Responses | 88 % Correct Responses | 92 % Correct Responses | 85 % Correct Responses |

URINE TOXICOLOGY REPORT SCREENING SECTION

Qualitative comparison of screening results available for each parameter.




- 1 Text section shows the correct response for the lab based on a comparison between the weighed in value and the lab's cut off value.
- 2 **Screening Results:** This chart is a quick visualisation of your performance over the last 20 samples. A result in the white section indicates a correct response. A result in the upper red section indicates a False Positive response, and a result in the lower red section indicates a False Negative response.
- 3 Comment section for RIQAS to provide your laboratory with additional relevant information regarding this sample, such as spiked metabolite concentration.
- 4 Screening result response categories. All abbreviations indicated at the bottom of the report page.
Key
TN - true negative TP - true positive FN - false negative
FP - false positive RC - sent for confirmation NT - not tested
- 5 **Screening Summary:** Your screening result shown in the appropriate response category and your cut off for this sample.
- 6 Screening results for all cut-offs returned for this sample within your method group.

- 7 Total screening results over all cut-offs for your laboratory's method.
- 8 Screening results for all cut-offs returned for this sample over all methods.
- 9 Total screening results over all cut-offs for all methods.
- 10 Screening results for other methods using same cut-off as your laboratory.
- 11 Performance history for this parameter, based on previous 10 samples.
- 12 Performance of your method over all cut-offs for this sample.
- 13 Performance history of your method over all cut-offs, based on the previous 10 samples.
- 14 Performance of all methods over all cut-offs for this sample.
- 15 Performance history of all methods over all cut-offs, based on the previous 10 samples.

URINE TOXICOLOGY REPORT QUANTITATIVE SECTION

Quantitative statistical comparison available for each parameter.



1 **Quantitative Text Section:** Comparison statistics. Caution is needed when the N value is too small to support statistical significance.

2 Your Result.

3 Your Mean for Comparison.

4 **Standard Deviation Index** = $\frac{\text{Your Result} - \text{Mean for Comparison}}{\text{SD of Mean for comparison}}$

5 Running mean SDI = average of last 10 SDIs for this parameter (If fewer than 10 results, "Too Few" is printed).

6 **Quantitative Results Histogram:** This graph provides a quick visualisation of how your quantitative result falls into the overall picture for all methods and your method group.

7 All available method statistics for this sample.

Your performance for each parameter is presented in a simple, convenient report.

Screening Results




- Categories are stated in your unit.
- Your method group.
- Your categories (available result options for chosen test strip and unit).
- All categories (result options) available for this parameter for any method (test strip).
- Results from all methods (test strips) returning results in the same categories as your lab.
- Results from all methods for all available categories.
- Your Result.
- Your Score:** Scores between 0-6 are acceptable, 7 borderline and 8 - 10 unacceptable.
- Target category and percentage of submitted results in that category.
- Performance Statement.

- Comments Box:** Provides number of correct scores and acceptable assessments for the last 6 samples.
- Categories Histogram:** A quick visualisation of how your lab's result falls into the overall picture for your categories.
- Possible reporting categories for your method.
- Your result is indicated by the black triangle.
- Levey-Jennings Chart:** Acceptable scores (0-6) have no shading, borderline scores (7) have light red shading, unacceptable scores (8-10) have dark red shading.
- Score for each sample number.
- Sample Number.
- Target Categories.
- All methods reported for this parameter.
- Detailed summary of results:** This table enables you to see how you compare to all other results.

SEROLOGY: SCREENING (QUALITATIVE) REPORT

Your performance for multiple samples is presented in a convenient single report per quarterly distribution.



1 Your qualitative result and chosen method are presented along with the acceptable result based on an 80% consensus. This consensus will be at the method level if there are >5 labs in the group or if there are <5 labs, will be at the all method level.

2 Overall Summary shows the number of results for this parameter and sample which are non-reactive, inconclusive or reactive.

3 Your Result is shown as a black triangle on the category chart compared to other laboratories in groups:

All Methods Your Method

4 Summary shows performance of all the methods used to analyse the parameter.

SEROLOGY: SCREENING (QUANTITATIVE) REPORT


Your performance for multiple samples is presented in a convenient single report per quarterly distribution.

Anti-Rubella IgG, IU/ml

Sample 2


| | N | Mean | CV% | U _m | SDPA | Exc. |
|------------------|-----|--------|------|----------------|-------|------|
| All methods | 210 | 92.574 | 37.2 | 2.97 | 34.42 | 31 |
| Abbott Architect | 39 | 83.219 | 8.7 | 1.46 | 7.27 | 5 |

| | | | |
|-----------------------|--------|-------|---------|
| ▲ Your Result | 84.800 | SDI | 0.22 |
| | | RMSDI | Too Few |
| ■ Mean for Comparison | 83.219 | | |



3

| Method | N | Mean | CV% | U _m |
|--------------------------------|----|---------|------|----------------|
| Biomerieux VIDAS | 48 | 150.979 | 9.8 | 2.97 |
| Abbott Architect | 44 | 83.219 | 8.7 | 1.46 |
| Roche Cobas 6000/8000 | 18 | 58.792 | 3.6 | 0.68 |
| Abbott Axsym | 17 | 108.206 | 18.0 | 6.09 |
| Siemens/DPC Immulite 2000/2500 | 17 | 90.800 | 6.2 | 1.94 |
| Roche Cobas 4000/e411 | 17 | 59.973 | 7.0 | 1.35 |
| Siemens/Bayer ADVIA Centaur | 14 | 120.775 | 11.0 | 5.88 |
| Roche Elecsys | 11 | 57.043 | 3.9 | 1.05 |
| Diasorin Liaison | 9 | 52.388 | 18.0 | 4.16 |
| Roche Modular E170 | 9 | 58.949 | 3.9 | 1.08 |
| Beckman DxI 600/800 | 6 | 125.817 | 7.4 | 4.75 |



4

1 Quantitative statistics for All Methods and Your Method are presented in your chosen unit along with your result and your performance scores (SDI and RMSDI).

2 Your result is presented on the bar graph as a black triangle, showing how you compare to:

All Methods Your Method

3 Multi Method Statistics section provides an easy way of assessing the performance of the methods used to analyse the parameter.

4 Levey-Jennings chart - Your SDIs for previous 20 samples.

QUANTITATIVE (END-OF-CYCLE REPORT)

The End-of-Cycle Report is sent to labs receiving standard reports at the end of each cycle and provides a complete summary of statistics. Results can also be compared to the previous cycle.

Albumin, g/l

Method: Bromocresol Purple
Instrument: Siemens/Dade Dimension RxL/Max/Xpand
Reagent: Siemens/Dade Behring

RIQAS TDPA: 7.1% **Biological Variation:** 3.9%

| Sample | Result | Unit | N | Mean for Comparison | CV% | Um | SDPA | SDI | TS | % Deviation |
|--------|--------|------|----|---------------------|-----|------|------|-------|-----|-------------|
| 1 | 28.200 | g/l | 68 | I 28.013 | 2.4 | 0.10 | 1.26 | 0.15 | 120 | 0.67 |
| 2 | 26.900 | g/l | 87 | I 26.853 | 2.7 | 0.10 | 1.21 | 0.04 | 120 | 0.17 |
| 3 | 39.900 | g/l | 71 | I 40.531 | 2.5 | 0.15 | 1.82 | -0.35 | 118 | -1.56 |
| 4 | 19.200 | g/l | 81 | I 19.429 | 2.5 | 0.07 | 0.87 | -0.26 | 120 | -1.18 |
| 5 | 41.700 | g/l | 67 | I 41.859 | 2.0 | 0.13 | 1.88 | -0.08 | 120 | -0.38 |
| 6 | 57.300 | g/l | 87 | I 57.257 | 2.7 | 0.21 | 2.58 | 0.02 | 120 | 0.08 |
| 7 | 45.000 | g/l | 72 | I 45.850 | 2.1 | 0.14 | 2.06 | -0.41 | 110 | -1.85 |
| 8 | 27.600 | g/l | 87 | I 28.013 | 2.5 | 0.09 | 1.26 | -0.33 | 120 | -1.47 |
| 9 | 41.200 | g/l | 70 | I 41.891 | 2.2 | 0.14 | 1.88 | -0.37 | 115 | -1.65 |
| 10 | 26.900 | g/l | 83 | I 26.742 | 3.3 | 0.12 | 1.20 | 0.13 | 120 | 0.59 |
| 11 | 40.700 | g/l | 71 | I 40.601 | 2.2 | 0.14 | 1.83 | 0.05 | 120 | 0.24 |
| 12 | 45.100 | g/l | 80 | I 45.456 | 2.2 | 0.14 | 2.04 | -0.17 | 120 | -0.78 |
| 13 | 27.300 | g/l | 63 | I 28.179 | 2.0 | 0.09 | 1.27 | -0.69 | 87 | -3.12 |

| | Cycle 45 | Cycle 46 |
|-----------------------------|----------|----------|
| Cycle Average SDI | -0.23 | -0.18 |
| Cycle Average TS | 110 | 116 |
| Cycle Average %DEV | -1.05 | -0.79 |
| Cycle Average Absolute SDI | 0.36 | 0.24 |
| Cycle Average Absolute %DEV | 1.63 | 1.06 |








CHART SECTION (END-OF-CYCLE REPORT)

Your results for current cycle shown in various diagrams.



- | | | |
|---|-----------------------------------|---|
| 1 | Levey-Jennings chart | <p>Shows your SDIs for a full cycle.</p> <ul style="list-style-type: none"> • Shows SDI (positive and negative) x Shows absolute SDI |
| 2 | Target Score chart | <p>Shows your Target Scores for a full cycle.</p> |
| 3 | %Deviation by sample chart | <p>Shows your %Deviations for a full cycle.</p> <p>Acceptable limits equal to TDPA unless alternative limits are registered by the lab.</p> <ul style="list-style-type: none"> • Shows %Deviation (positive and negative) x Shows absolute %Deviation |
| 4 | %Deviation by Concentration chart | <p>Shows your results for a full cycle.</p> |

TEXT SECTION (END-OF-CYCLE REPORT)

The text section summarises the statistical information for all samples.

1 Albumin, g/l

2 **Method:** Bromocresol Purple
Instrument: Siemens/Dade Dimension RxL/Max/Xpand
Reagent: Siemens/Dade Behring

3 **RIQAS TDPA:** 7.1% **Biological Variation:** 3.9%

Your assay details at the end of the cycle. The RIQAS TDPA and biological variation for the parameter are shown if available.

4 5 6 7 8 9 10 11 12 13 14

| Sample | Result | Unit | N | Mean | SDPA | U _m | CV% | SDI | TS | % Deviation |
|--------|--------|------|----|----------|------|----------------|-----|-------|-----|-------------|
| 1 | 28.200 | g/l | 68 | I 28.013 | 1.26 | 0.10 | 2.4 | 0.15 | 120 | 0.7 |
| 2 | 26.900 | g/l | 87 | I 26.853 | 1.21 | 0.10 | 2.7 | 0.04 | 120 | 0.2 |
| 3 | 39.900 | g/l | 71 | M 40.531 | 1.82 | 0.15 | 2.5 | -0.36 | 116 | -1.5 |
| 4 | 19.200 | g/l | 81 | I 19.429 | 0.87 | 0.07 | 2.5 | -0.27 | 120 | -1.2 |
| 5 | 41.700 | g/l | 67 | I 41.942 | 1.88 | 0.13 | 2.0 | -0.09 | 120 | -0.4 |
| 6 | 57.300 | g/l | 87 | I 57.257 | 2.58 | 0.21 | 2.7 | 0.02 | 120 | 0.1 |
| 7 | 45.000 | g/l | 72 | I 45.850 | 2.06 | 0.14 | 2.1 | -0.43 | 108 | -1.8 |
| 8 | 27.600 | g/l | 87 | I 28.011 | 1.26 | 0.09 | 2.5 | -0.34 | 118 | -1.5 |
| 9 | 41.200 | g/l | 70 | I 41.823 | 1.88 | 0.14 | 2.2 | -0.38 | 113 | -1.6 |
| 10 | 26.900 | g/l | 83 | I 26.742 | 1.20 | 0.12 | 3.3 | 0.14 | 120 | 0.6 |
| 11 | 40.700 | g/l | 71 | I 40.601 | 1.83 | 0.13 | 2.2 | 0.06 | 120 | 0.2 |
| 12 | 45.100 | g/l | 80 | I 45.119 | 2.05 | 0.14 | 2.2 | -0.18 | 120 | -0.8 |
| 13 | 27.300 | g/l | 63 | I 28.454 | 1.27 | 0.09 | 2.0 | -0.72 | 86 | -3.1 |

Summary of your results and statistics are shown, including Mean for Comparison, SDPA, %CV, U_m, SDI, Target Score, %Deviation.

| | Cycle 45 | Cycle 46 |
|--------------------------------------|----------|----------|
| 15 Cycle Average SDI | -0.23 | -0.18 |
| Cycle Average TS | 110 | 116 |
| Cycle Average %DEV | -1.05 | -0.79 |
| 16 Cycle Average Absolute SDI | 0.36 | 0.24 |
| Cycle Average Absolute %DEV | 1.63 | 1.06 |

Table containing a summary of your performance for previous cycle and current cycle, including Average Absolute SDIs and %Deviations.

TEXT SECTION (END-OF-CYCLE REPORT)

- 1 Report presented in your chosen unit
- 2 Your assay details as of the last sample
- 3 RIQAS TDPA and Biological variation
- 4 Sample number
- 5 Your results for each sample
- 6 Unit your result was returned in
- 7 Number of results used for statistical analysis
- 8 Mean for Comparison (including comparison level)
- 9 SDPA = Standard Deviation for performance assessment
- 10 Uncertainty of Mean for Comparison
- 11 Coefficient of Variation (%)
- 12 Your Standard Deviation Index
- 13 Your Target Score
- 14 Your %Deviation

- 15 Cycle average of your performance indicators – Standard Deviation Index, Target Score and %Deviation.

$$\text{Cycle Average SDI} = \frac{\text{(Sum of SDIs returned for the completed cycle)}}{\text{(Number of samples returned in cycle)}}$$

$$\text{Cycle Average Target Score} = \frac{\text{(Sum of your Target Scores returned for the completed cycle)}}{\text{(Number of samples returned in cycle)}}$$

$$\text{Cycle Average \%Deviation} = \frac{\text{(Sum of your \%Deviations returned for the completed cycle)}}{\text{(Number of samples returned in cycle)}}$$


- 16 Cycle average for Absolute values of your SDI and %Deviation. Absolute values show how far a value is from zero regardless of the sign. This is an indication of the magnitude of accuracy.

$$\text{Cycle Average Absolute SDI} = \frac{\text{(Sum of your Absolute SDIs returned for the completed cycle)}}{\text{(Number of samples returned in cycle)}}$$

$$\text{Cycle Average Absolute \%Deviation} = \frac{\text{(Sum of your Absolute \%Deviations returned for the completed cycle)}}{\text{(Number of samples returned in cycle)}}$$

CURRENT & PREVIOUS CYCLE ABSOLUTE SDIs (END-OF-CYCLE REPORT)

Based on the cycle average absolute SDI, this chart provides a visual representation of your laboratory's performance compared to the previous cycle.



- 1** Parameter list

List of all parameters registered.
- 2** Results for previous cycle

Indicated by open circle on the chart.
- 3** Report title - Cycle Average Absolute SDI

This shows your performance this cycle compared to the previous cycle.
- 4** Results for current cycle

Indicated by a closed circle on the chart.
- 5** Legend

Cycle Average Absolute SDIs are shown for:

 - L** Your results throughout the cycle
 - C** All labs within your own country
 - W** All labs Worldwide
- 6** Graphical representation of Absolute SDIs

Acceptable performance is ≤ 2 .

If Absolute SDI for current cycle is less than that for the previous cycle, this is indicated by a green circle.

If Absolute SDI for current cycle is greater than that for the previous cycle, this is indicated by a red circle.

The closer the circle is to zero, the better the performance.

CERTIFICATE OF PERFORMANCE (END-OF-CYCLE REPORT)

An End-of-Cycle report will be issued for all registrations. However, the Certificate of Performance will only be available for parameters where results for at least 50% of samples in the cycle have been returned. Labs joining after the beginning of the cycle will only receive the Certificate of Performance if they meet this criterion. Any parameters not included on the Certificate of Acceptable Performance will be listed on the Notification of Unacceptable Performance.

RIQAS  **RANDOX INTERNATIONAL QUALITY ASSESSMENT SCHEME**

CERTIFICATE OF ACCEPTABLE PERFORMANCE

RIQAS Department
Randox Laboratories
CRUMLIN
COUNTY ANTRIM
BT29 4QY
UNITED KINGDOM

1

2 — LABORATORY REF. NO. XX/X

3 — CLINICAL CHEMISTRY - CYCLE 47

4 — 11/03/2013

This is to certify that the above participant took part in a cycle of external quality assessment and achieved an acceptable level of performance (Cycle Average Absolute SDI ≤ 2) for the following parameters:

5

6 — Cycle Average Absolute SDI

| | |
|--|------|
| Albumin - Bromocresol Purple - Siemens/Dade Dimension RxL/Max/Xpand | 0.50 |
| Alkaline Phosphatase - Dade Dimension, AMP buffer - Siemens/Dade Dimension RxL/Max/Xpand | 1.22 |
| ALT (GPT) - Tris buffer with P5P - Siemens/Dade Dimension RxL/Max/Xpand | 0.53 |
| Amylase, Total - Dade Behring 2-chloro-pNPG3 - Siemens/Dade Dimension RxL/Max/Xpand | 0.34 |
| AST (GOT) - Tris buffer with P5P - Siemens/Dade Dimension RxL/Max/Xpand | 0.55 |
| Bicarbonate - Enzymatic - Siemens/Dade Dimension RxL/Max/Xpand | 1.08 |
| Bilirubin, Direct - Diazo with Sulphanilic Acid - Siemens/Dade Dimension RxL/Max/Xpand | 0.19 |
| Bilirubin, Total - Diazo with Sulphanilic Acid - Siemens/Dade Dimension RxL/Max/Xpand | 0.26 |
| Calcium - Cresolphthalein complexone - Siemens/Dade Dimension RxL/Max/Xpand | 0.49 |
| Chloride - ISE, indirect - Siemens/Dade Dimension RxL/Max/Xpand | 0.70 |
| Cholesterol - Dimension-Dade Behring reagents - Siemens/Dade Dimension RxL/Max/Xpand | 0.54 |
| CK, Total - CK-NAC (IFCC) - Siemens/Dade Dimension RxL/Max/Xpand | 0.26 |
| Creatinine - Alkaline picrate no deprot. - Siemens/Dade Dimension RxL/Max/Xpand | 0.44 |
| GGT - Gamma glut'3-carb'4-nitro (IFCC) - Siemens/Dade Dimension RxL/Max/Xpand | 0.25 |
| Glucose - Hexokinase - Siemens/Dade Dimension RxL/Max/Xpand | 0.70 |

| | | |
|---|---------------------------|---|
| 1 | Full registration address | Your full registration address details. |
| 2 | Your lab reference number | Used to identify each lab. |
| 3 | Programme / cycle number | Programme and current, completed cycle number. |
| 4 | Date | Date End-of-Cycle report is issued. |
| 5 | Parameters | List of parameters including the assay details for which cycle absolute SDI is ≤ 2 . |
| 6 | Average Absolute SDI | Your Cycle Average Absolute SDI. |

MONITORING EQA PERFORMANCE


Each EQA report should be evaluated and any poor performance investigated. A step by step approach should be adopted consisting of the following three steps:

1. Investigate the source of the problem

In order to identify the source of the problem, it is useful to be aware of the most common causes of poor EQA performance. Errors can occur at any stage of the testing process; however, EQA is most concerned with detecting analytical errors i.e. errors that occur during the analysis of the sample.

Most analytical errors can be easily divided into three main areas; clerical errors, systematic errors and random errors. Systematic errors result in inaccurate results that consistently show a positive or negative bias. Random errors, on the other hand, affect precision and result in fluctuations in either direction.

It may be possible that, after extensive investigations, the root cause of the poor performance cannot be established. Poor performance for a single sample could be attributed to random error. If poor performance has been noted for several samples, a systematic error is the most likely cause and the analytical process should be reviewed.



Clerical errors

- Transcription errors
- Incorrect units used
- Incorrect sample tested
- Incorrect method classification
- Calculation/conversion error

Systematic errors

- Sample/Reagent prep/handling
- Reagent/calibrator/standardisation change
- Instrument/reagent/calibrator fault
- Inexperienced operators
- Reagent deterioration
- Inappropriate method

Random errors

- Bubbles in reagent
- Bubbles in reagent/sample pipette
- Temperature fluctuations
- Poor pipetting technique
- Poor operator technique

The flowchart (page 29) is designed to help you investigate any apparent poor performance.

2. Implement corrective actions

Some errors can be readily recognised as simple clerical errors and easily corrected. If there is evidence of systematic or random error however more detailed corrective actions must be taken.

Systematic Error

In the event of a systematic error, the following suggested actions may help to resolve the problem:

- Perform instrument maintenance
- Recalibrate instrument
- Review reagent/sample storage
- Check pipettes
- Prepare fresh reagents & re-run sample
- Perform staff training

Random Error

If all possible causes have been excluded, a single unacceptable result is most likely due to random error. Re-run the sample; if the result of repeat analysis is acceptable then corrective action is not required. If the issue persists, investigate possible sources of systematic error.

3. Check the effectiveness of corrective actions

The effectiveness or impact of any corrective actions taken can be assessed by continuing to monitor analytical performance over time.

MONITORING EQA PERFORMANCE

A checklist similar to the one below is extremely useful when investigating poor EQA performance and may help you to determine the root cause of the problem and initiate corrective actions.

Laboratory:
 Cycle Number: Sample Number:
 Analysis Date: Analyte:
 Mean for Comparison: Lab Result: SDI: %Dev:

1. Specimen Handling

- a. Samples received in good condition Y N
- b. Samples stored/prepared appropriately Y N
- c. Integrity of the sample is acceptable Y N

2. Clerical

- a. Correct result entered Y N
- b. Correct use of decimal point and units Y N
- c. Calculations, if any, performed correctly (even if automated) Y N
- d. Conversion factors applied to results before submission Y N

3. Registration and Mean for Comparison

- a. Registered in the correct method/instrument group Y N
- b. Changed method or instrument without advising RIQAS Y N
- c. Peer Group changed due to the number of participants returning results e.g. from method to instrument Y N
- d. An obvious bias between method and instrument means (check histogram and stats sections) Y N

4. Internal Quality Control

- a. %Deviation of IQC (at similar conc to that of EQA) on sample analysis date acceptable Y N
- b. Shift in IQC in the periods just before and after EQA sample analysis Y N
- c. Trends in IQC in the periods before and after EQA sample analysis Y N
- d. Random IQC variation on sample analysis date Y N

- e. Error due to imprecision; check IQC in terms of %Deviation compared to deviation observed in EQA Y N
- f. IQC target correctly assigned Y N

5. Calibration

- a. Date of last calibration
- b. Calibration frequency acceptable Y N
- c. Last calibration acceptable Y N

6. Instrument

- a. Daily maintenance performed on date of sample analysis Y N
- b. Special maintenance performed prior to sample analysis Y N
- c. Instrument operated correctly Y N
- d. Operator fully trained Y N

7. Reagents

- a. Reagents prepared and stored correctly Y N
- b. Reagents within open vial stability Y N

8. EQA sample

- a. Initial value
- b. Re-run value
- c. Issue observed in previous EQA samples at a similar concentration (check %Deviation by concentration and Levey Jennings charts) Y N
- d. All parameters affected (to the same extent) - possible reconstitution error (check %Deviation on summary pages) Y N

Conclusion:


Remedial Action:

Lab Manager: Date:

Lab Director: Date:

MONITORING EQA PERFORMANCE

The flow chart below can be used to help identify a possible root cause for poor EQA performance.



Ammonia/Ethanol Programme+ *With target scoring*

RQ9164 (2 ml)
2 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

Ammonia Ethanol

Anti-TSH Receptor Programme+ *With target scoring*

RQ9174 (1 ml)
1 Parameter
Samples every month, 1 x 12 month cycle, 12 month subscription

Anti-TSH Receptor (TRAb)

Blood Gas Programme *With target scoring*

| | |
|---|--|
| RQ9134 (1.8 ml) First registered instrument 11 Parameters Samples every month, 1 x 12 month cycle, 12 month subscription | RQ9134/A (1.8 ml) Subsequent instruments 11 Parameters |
|---|--|

| | | | |
|------------------|-------------------------|------------------|-----------------|
| Bicarbonate | CO ₂ (Total) | K+ | pH |
| Ca ⁺⁺ | Glucose | Na+ | pO ₂ |
| Cl- | Lactate | pCO ₂ | |

BNP Programme+ *With target scoring*

RQ9165 (1 ml)
1 Parameter
Samples every month, 1 x 12 month cycle, 12 month subscription

BNP

Cardiac Programme *With target scoring*

| | | |
|--|--------------------------------------|---|
| RQ9127/a (1 ml) 2 Parameters only (choose from 7) Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription | RQ9127/b (1 ml) Full 7 Parameters | RQ9186 (1ml) Full 7 Parameters Samples every month, 1 x 12 monthly cycle, 12 month subscription |
|--|--------------------------------------|---|

| | | | |
|-------------------------------|------------------------------|-------------------------|------------|
| CK, Total CK-MB (Activity) | CK-MB (Mass) Homocysteine | Myoglobin Troponin I | Troponin T |
|-------------------------------|------------------------------|-------------------------|------------|

Cardiac Plus Programme • **coming in 2021*

RQ9190 (3 ml)
11 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

| | | | |
|---|------------------------------------|---------------------------------|--------------------------|
| CK, Total CK-MB Activity CK-MB Mass | D-dimer Digoxin Homocysteine | hsCRP Myoglobin NT proBNP | Troponin I Troponin T |
|---|------------------------------------|---------------------------------|--------------------------|

Cerebrospinal Fluid Programme+ *With target scoring*

RQ9168 (3 ml)
7 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

| | | | |
|---------------------|----------------|----------------------------|--------|
| Albumin Chloride | Glucose IgG | Lactate Protein (Total) | Sodium |
|---------------------|----------------|----------------------------|--------|

RIQAS PROGRAMMES

Coagulation Programme *With target scoring*

RQ9135/a (1 ml)
5 Selected parameters only + 1 pilot
(aPTT, PT, TT, Fibrinogen, Antithrombin III)
Samples every month, 1 x 12 month cycle, 12 month subscription

aPTT
PT (including INR)
TT
Fibrinogen
Antithrombin III

RQ9135/b (1 ml)
Full 16 Parameters + 1 pilot

D-dimer*
Factor II
Factor V
Factor VII
Factor VIII

Factor IX
Factor X
Factor XI
Factor XII
Plasminogen

Protein C
Protein S

CO-Oximetry Programme+

RQ9177 (1.2 ml)
First registered instrument
7 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

Carboxyhaemoglobin (COHb / HbCO)
Deoxyhaemoglobin (HHb)

Methaemoglobin (MetHb)
Oxygen Content (O2CT)

Oxygen Saturation (sO2 / Vol O2)
Oxyhaemoglobin (O2Hb / HbO2)

Total Haemoglobin (tHb)

CYFRA 21-I Programme+

RQ9175 (1 ml)
1 Parameter
Samples every month, 1 x 12 month cycle, 12 month subscription

CYFRA 21-I (Cytokeratin 19)

ESR Programme+

RQ9163 (4.5 ml)
1 Parameter
2 samples per quarterly distribution, 1 x 12 month cycle, 12 month subscription

ESR (Erythrocyte Sedimentation Rate)

General Clinical Chemistry Programme *With target scoring*

RQ9112/a (5 ml)
10 Parameters + 4 pilots

RQ9112/b (5 ml)
17 Parameters + 4 pilots

RQ9112/c (5 ml)
Full 52 Parameters + 4 pilots

RQ9128 (5ml)
Full 52 Parameters + 4 pilots
Samples every month, 1 x 12 monthly cycle, 12 month subscription

Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription, reference method values

ACE (Angiotensin Converting Enzyme)
Acid Phosphatase (Prostatic)
Acid Phosphatase (Total)
Albumin
Alkaline Phosphatase
ALT (ALAT)
Amylase (Pancreatic)
Amylase (Total)
AST (ASAT)
Bicarbonate
Bile Acids
Bilirubin (Direct)
Bilirubin (Total)
Calcium
Calcium, Adjusted*

Calcium (Ionised)
Chloride
Cholesterol
Cholinesterase
CK, Total (CPK)
Copper
Creatinine
D-3-Hydroxybutyrate
eGFR (estimated glomerular filtration rate)*
Fructosamine
γGT
GLDH
Glucose
HBDH
HDL-Cholesterol

Iron
Lactate
LD (LDH)
LDL-Cholesterol*
Lipase
Lithium
Magnesium
NEFA
Non-HDL Cholesterol*
Osmolality
Phosphate (Inorganic)
Potassium
Protein (Total)
PSA
Sodium

TIBC
T₃ (Free)
T₃ (Total)
T₄ (Free)
T₄ (Total)
Triglycerides
TSH
UIBC
Urea
Uric Acid
Zinc

Glycated Haemoglobin Programme (HbA1c) *With target scoring*

RQ9129 (0.5ml)
2 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

HbA1c

Total Haemoglobin

 = Liquid ready-to-use samples  = Lyophilised samples **PURPLE** = The only parameters available on RQ9135/a + = Not accredited * = Pilot study ongoing • = Accreditation status pending

Haematology Programme *With target scoring*

RQ9118 (2 ml)

11 Parameters

Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription

Haematocrit (HCT)
Haemoglobin (Hb)
Mean Cell Haemoglobin (MCH)

Mean Cell Haemoglobin Concentration (MCHC)
Mean Cell Volume (MCV)
Mean Platelet Volume (MPV)

RQ9140 (2ml)

11 Parameters

Samples every month, 1 x 12 monthly cycle, 12 month subscription

Platelets (PLT)
Plateletcrit (PCT)
Red Blood Cell Count (RBC)

Red Cell Distribution Width (RDW)
Total White Blood Cell Count (WBC)

Human Urine Programme *With target scoring*

RQ9115 (10 ml)

25 Parameters

Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription

ACR
Albumin/Microalbumin
Amylase
Calcium
Chloride
Copper
Cortisol

Creatinine
Dopamine
Epinephrine
Glucose
Metanephrine
Norepinephrine

RQ9185 (10ml)

25 Parameters

Samples every month, 1 x 12 monthly cycle, 12 month subscription

Normetanephrine
Magnesium
Osmolality
Oxalate
Phosphate (Inorganic)
Potassium

Protein (Total)
Sodium
Urea
Uric Acid
VMA
5-HIAA

Immunoassay Programme *With target scoring*

RQ9125/a (5 ml)

4 Parameters only + 2 pilots

Samples every two weeks, 2 x 6 monthly cycles, 12 month subscription (RQ9125/a, RQ9125/b, RQ9125/c)

ACTH
AFP
Aldosterone
Amikacin
Androstenedione
 β -2-Microglobulin
CA125
CA15-3
CA19-9
Carbamazepine
CEA
Cortisol
C-Peptide

RQ9125/b (5 ml)

13 Parameters only + 2 pilots

DHEA-Sulphate
DHEA Unconjugated
Digoxin
Ferritin
Folate
FSH
Gentamicin
GH
hCG
IgE
Insulin
LH
Oestradiol

RQ9125/c (5 ml)

Full 49 Parameters + 2 pilots

17-OH-Progesterone
Paracetamol
Phenobarbital
Phenytoin
Progesterone
Prolactin
PSA (Free)
PSA (Total)
PTH
Salicylate
SHBG
 T_3 (Free)
 T_3 (Total)

RQ9130 (5 ml)

Full 49 Parameters + 2 pilots
Samples every month, 1 x 12 month cycle, 12 month subscription (RQ9130)

T_4 (Free)
 T_4 (Total)
Testosterone (Free)*
Testosterone (Total)
Theophylline
Thyroglobulin
TSH
Valproic Acid
Vancomycin
Vitamin B12
1-25-(OH) $_2$ -Vitamin D*
25-OH-Vitamin D

Immunoassay Speciality 1 Programme+ *With target scoring*

RQ9141 (2 ml)

9 Parameters + 1 pilot

Samples every month, 1 x 12 month cycle, 12 month subscription

1-25-(OH) $_2$ -Vitamin D*
25-OH-Vitamin D
C-Peptide

Anti-TG
Anti-TPO
IGF-I

Osteocalcin
Procalcitonin
PTH

Insulin

Immunoassay Speciality 2 Programme+ *With target scoring*

RQ9142 (1 ml)

5 Parameters

Samples every month, 1 x 12 month cycle, 12 month subscription

Calcitonin
Gastrin

Procalcitonin

Plasma Renin Activity

Renin (Direct Concentration)

Immunosuppressant Programme+

RQ9159 (2 ml)

4 Parameters

Samples every month, 1 x 12 month cycle, 12 month subscription, reference method values

Ciclosporin

Everolimus

Sirolimus

Tacrolimus

 = Liquid ready-to-use samples

 = Lyophilised samples

PURPLE = The only parameters available on RQ9135/a

+ = Not accredited

* = Pilot study ongoing

• = Accreditation status pending

RIQAS PROGRAMMES

Lipid Programme *With target scoring*

RQ9126/a (3 ml)
3 Parameters only (choose from 7)
Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription

Apolipoprotein A I
Apolipoprotein B

RQ9126/b (3 ml)
Full 7 Parameters

Cholesterol (Total)
HDL-Cholesterol

LDL-Cholesterol
Lipoprotein (a)

Triglycerides

Maternal Screening Programme *With target scoring*

RQ9137 (1 ml)
6 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

AFP
free β -hCG

Total hCG
Inhibin A

PAPP-A

Unconjugated Oestriol

Serology (EBV) Programme+

RQ9153 (1 ml)
3 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription, Quantitative and Qualitative results

Anti-EBV VCA IgG

Anti-EBNA IgG

Anti-EBV VCA IgM

Serology (HIV-Hepatitis) Programme+

RQ9151 (1.8 ml)
10 Parameters + 6 pilots
Samples every month, 1 x 12 month cycle, 12 month subscription, Quantitative and Qualitative results

Anti-CMV (Total)
Anti-HAV IgM*
Anti-HAV (Total)*
Anti-HBc

Anti-HBc IgM*
Anti-HBe (Total)*
Anti-HBs (Total)*
Anti-HCV

Anti-HIV-1
Anti-HIV-2
Anti-HIV combined
Anti-HTLV I

Anti-HTLV II
Anti-HTLV combined
HBsAg
P24*

Serology (Syphilis) Programme+

RQ9154 (1 ml)
1 Parameter
Samples every month, 1 x 12 month cycle, 12 month subscription, Quantitative and Qualitative results

Syphilis (Methods available include immunoassay RPR, VDRL and TPHA)

Serology (ToRCH) Programme+

RQ9152 (1 ml)
12 Parameters + 3 pilots
Samples every month, 1 x 12 month cycle, 12 month subscription, Quantitative and Qualitative results

Anti-CMV IgG
Anti-CMV IgM
Anti-HSV1 IgG
Anti-HSV1 IgM

Anti-HSV2 IgG
Anti-HSV2 IgM
Anti-HSV1/2 IgG
Anti-HSV1/2 IgM

Anti-Measles IgG*
Anti-Mumps IgG*
Anti-Rubella IgG
Anti-Rubella IgM

Anti-Toxoplasma IgG
Anti-Toxoplasma IgM
Anti-VZV IgG*

Specific Proteins Programme *With target scoring*

RQ9114 (3 ml)
26 Parameters
Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription

AFP
Albumin
 α -1-Acid glycoprotein
 α -1-Antitrypsin
 α -2-Macroglobulin
Anti Streptolysin O
Antithrombin III

β -2-Microglobulin
Ceruloplasmin
Complement C₃
Complement C₄
C-Reactive Protein
Ferritin
Haptoglobin

RQ9187 (1ml)
26 Parameters
Samples every month, 1 x 12 monthly cycle, 12 month subscription

IgA
IgE
IgG
IgM
Kappa Light Chain (Free)
Kappa Light Chain (Total)
Lambda Light Chain (Free)

Lambda Light Chain (Total)
Prealbumin (Transthyretin)
Retinol Binding Protein
Rheumatoid Factor
Transferrin

 = Liquid ready-to-use samples  = Lyophilised samples **PURPLE** = The only parameters available on RQ9135/a + = Not accredited * = Pilot study ongoing • = Accreditation status pending

Sweat Testing Programme+

RQ9173 (2 ml)
2 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

Chloride Conductivity

Therapeutic Drugs Programme *With target scoring*

RQ9111 (5 ml)
18 Parameters
Samples every 2 weeks, 2 x 6 monthly cycles, 12 month subscription, Weighed-in values

| | | | |
|---------------|-----------------------------|----------------|---------------|
| Amikacin | Ethosuximide | Phenobarbital | Tobramycin |
| Caffeine | Gentamicin | Phenytoin | Valproic Acid |
| Carbamazepine | Lithium | Primidone | Vancomycin |
| Ciclosporin | Methotrexate | Salicylic Acid | |
| Digoxin | Paracetamol (Acetaminophen) | Theophylline | |

Trace Elements In Blood Programme+

RQ9172 (3 ml)
7 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

Copper Lead Manganese Zinc
Iodine Magnesium Selenium

Trace Elements In Serum Programme+

RQ9170 (3 ml)
10 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

Aluminium Copper Manganese Zinc
Chromium Iodine Nickel
Cobalt Lead Selenium

Trace Elements In Urine Programme+

RQ9171 (3 ml)
11 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

Cadmium Copper Magnesium Nickel
Chromium Iodine Manganese Thallium
Cobalt Lead Molybdenum

Urinalysis Programme+ *With scoring*

RQ9138 (12 ml)
14 Parameters
Samples every 2 months, 1 x 12 month cycle, 12 month subscription

Albumin Galactose Leucocytes Specific Gravity
Bilirubin Glucose Nitrite Urobilinogen
Blood hCG pH
Creatinine Ketones Protein

Urine Toxicology Programme+

RQ9139 (5 ml)
20 Parameters
Samples every month, 1 x 12 month cycle, 12 month subscription

| | | | |
|--------------------|-------------------|-----------------|---------------|
| Benzoylgonine | d-Methamphetamine | MDMA | Phenobarbital |
| Buprenorphine | EDDP | Methadone | Secobarbital |
| Cannabinoids (THC) | Ethanol | Nortriptyline | |
| Cotinine | Free Morphine | Norpropoxyphene | |
| Creatinine | Lorazepam | Oxazepam | |
| d-Amphetamine | LSD | Phencyclidine | |

 = Liquid ready-to-use samples  = Lyophilised samples **PURPLE** = The only parameters available on RQ9135/a + = Not accredited * = Pilot study ongoing • = Accreditation status pending

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| # | Parameter | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-I + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality I + |
|---|-------------------------------------|---------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|
| | 1-25-(OH) ₂ -Vitamin D* | | | | | | | | | | | | | | | | X | X |
| | 17-OH-Progesterone | | | | | | | | | | | | | | | | X | |
| | 25-OH-Vitamin D | | | | | | | | | | | | | | | | X | X |
| | 5-HIAA | | | | | | | | | | | | | | | X | | |
| A | α-1-Acid Glycoprotein | | | | | | | | | | | | | | | | | |
| | α-1-Antitrypsin | | | | | | | | | | | | | | | | | |
| | α-2-Macroglobulin | | | | | | | | | | | | | | | | | |
| | ACE (Angiotensin Converting Enzyme) | | | | | | | | | | | | X | | | | | |
| | Acid Phosphatase (Prostatic) | | | | | | | | | | | | X | | | | | |
| | Acid Phosphatase (Total) | | | | | | | | | | | | X | | | | | |
| | ACR | | | | | | | | | | | | | | | X | | |
| | ACTH | | | | | | | | | | | | | | | | | X |
| | AFP | | | | | | | | | | | | | | | | | X |
| | Albumin | | | | | | | X | | | | | X | | | X | | |
| | Aldosterone | | | | | | | | | | | | | | | | | X |
| | Alkaline Phosphatase | | | | | | | | | | | | X | | | | | |
| | ALT (ALAT) | | | | | | | | | | | | X | | | | | |
| | Aluminium | | | | | | | | | | | | | | | | | |
| | Amikacin | | | | | | | | | | | | | | | | | X |
| | Ammonia | | X | | | | | | | | | | | | | | | |
| | Amylase (Pancreatic) | | | | | | | | | | | | | X | | | | |
| | Amylase (Total) | | | | | | | | | | | | | X | | X | | |
| | Androstenedione | | | | | | | | | | | | | | | | | X |
| | Anti Streptolysin O (ASO) | | | | | | | | | | | | | | | | | |
| | Anti-CMV | | | | | | | | | | | | | | | | | |
| | Anti-CMV IgG | | | | | | | | | | | | | | | | | |
| | Anti-CMV IgM | | | | | | | | | | | | | | | | | |
| | Anti-EBNA IgG | | | | | | | | | | | | | | | | | |
| | Anti-EBV VCA IgG | | | | | | | | | | | | | | | | | |
| | Anti-EBV VCA IgM | | | | | | | | | | | | | | | | | |
| | Anti-HAV IgM* | | | | | | | | | | | | | | | | | |
| | Anti-HAV (Total)* | | | | | | | | | | | | | | | | | |
| | Anti-HBc | | | | | | | | | | | | | | | | | |
| | Anti-HBc IgM* | | | | | | | | | | | | | | | | | |
| | Anti-HBe (Total)* | | | | | | | | | | | | | | | | | |
| | Anti-HBs (Total)* | | | | | | | | | | | | | | | | | |
| | Anti-HCV | | | | | | | | | | | | | | | | | |
| | Anti-HIV-1 | | | | | | | | | | | | | | | | | |
| | Anti-HIV-1 & 2 Combined | | | | | | | | | | | | | | | | | |
| | Anti-HIV-2 | | | | | | | | | | | | | | | | | |
| | Anti-HSV-1 & 2 IgG Combined | | | | | | | | | | | | | | | | | |

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunoassay Speciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | | |
|----------------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|--|-------------------------------------|---|
| | | | | | | | | | | | | | | | | | I-25-(OH) ₂ -Vitamin D* | # |
| | | | | | | | | | | | | | | | | | 17-OH-Progesterone | |
| | | | | | | | | | | | | | | | | | 25-OH-Vitamin D | |
| | | | | | | | | | | | | | | | | | 5-HIAA | |
| | | | | | | | | X | | | | | | | | | α-1-Acid Glycoprotein | A |
| | | | | | | | | X | | | | | | | | | α-1-Antitrypsin | |
| | | | | | | | | X | | | | | | | | | α-2-Macroglobulin | |
| | | | | | | | | | | | | | | | | | ACE (Angiotensin Converting Enzyme) | |
| | | | | | | | | | | | | | | | | | Acid Phosphatase (Prostatic) | |
| | | | | | | | | | | | | | | | | | Acid Phosphatase (Total) | |
| | | | | | | | | | | | | | | | | | ACR | |
| | | | | | | | | | | | | | | | | | ACTH | |
| | | | X | | | | | X | | | | | | | | | AFP | |
| | | | | | | | | X | | | | | | | X | | Albumin | |
| | | | | | | | | | | | | | | | | | Aldosterone | |
| | | | | | | | | | | | | | | | | | Alkaline Phosphatase | |
| | | | | | | | | | | | | | | | | | ALT (ALAT) | |
| | | | | | | | | | | | | X | | | | | Aluminium | |
| | | | | | | | | | | X | | | | | | | Amikacin | |
| | | | | | | | | | | | | | | | | | Ammonia | |
| | | | | | | | | | | | | | | | | | Amylase (Pancreatic) | |
| | | | | | | | | | | | | | | | | | Amylase (Total) | |
| | | | | | | | | | | | | | | | | | Androstenedione | |
| | | | | | | | | X | | | | | | | | | Anti Streptolysin O (ASO) | |
| | | | | | X | | | | | | | | | | | | Anti-CMV | |
| | | | | | | | X | | | | | | | | | | Anti-CMV IgG | |
| | | | | | | | X | | | | | | | | | | Anti-CMV IgM | |
| | | | | X | | | | | | | | | | | | | Anti-EBNA IgG | |
| | | | | X | | | | | | | | | | | | | Anti-EBV VCA IgG | |
| | | | | X | | | | | | | | | | | | | Anti-EBV VCA IgM | |
| | | | | | X | | | | | | | | | | | | Anti-HAV IgM* | |
| | | | | | X | | | | | | | | | | | | Anti-HAV (Total)* | |
| | | | | | X | | | | | | | | | | | | Anti-HBc | |
| | | | | | X | | | | | | | | | | | | Anti-HBc IgM* | |
| | | | | | X | | | | | | | | | | | | Anti-HBe (Total)* | |
| | | | | | X | | | | | | | | | | | | Anti-HBs (Total)* | |
| | | | | | X | | | | | | | | | | | | Anti-HCV | |
| | | | | | X | | | | | | | | | | | | Anti-HIV-1 | |
| | | | | | X | | | | | | | | | | | | Anti-HIV-1 & 2 Combined | |
| | | | | | X | | | | | | | | | | | | Anti-HIV-2 | |
| | | | | | | X | | | | | | | | | | | Anti-HSV-1 & 2 IgG Combined | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| | | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-I + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality + |
|---|-----------------------------|---------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|
| A | Anti-HSV-1 & 2 IgM Combined | | | | | | | | | | | | | | | | | |
| | Anti-HSV1 IgG | | | | | | | | | | | | | | | | | |
| | Anti-HSV1 IgM | | | | | | | | | | | | | | | | | |
| | Anti-HSV2 IgG | | | | | | | | | | | | | | | | | |
| | Anti-HSV2 IgM | | | | | | | | | | | | | | | | | |
| | Anti-HTLV-I & 2 Combined | | | | | | | | | | | | | | | | | |
| | Anti-HTLV-I | | | | | | | | | | | | | | | | | |
| | Anti-HTLV-II | | | | | | | | | | | | | | | | | |
| | Anti-Measles IgG* | | | | | | | | | | | | | | | | | |
| | Anti-Mumps IgG* | | | | | | | | | | | | | | | | | |
| | Anti-Rubella IgG | | | | | | | | | | | | | | | | | |
| | Anti-Rubella IgM | | | | | | | | | | | | | | | | | |
| | Anti-TG | | | | | | | | | | | | | | | | | X |
| | Anti-VZV IgG* | | | | | | | | | | | | | | | | | |
| | Antithrombin III | | | | | | | | | X | | | | | | | | |
| | Anti-Toxoplasma IgG | | | | | | | | | | | | | | | | | |
| | Anti-Toxoplasma IgM | | | | | | | | | | | | | | | | | |
| | Anti-TPO | | | | | | | | | | | | | | | | | X |
| | Anti-TSH Receptor (TRAb) | | X | | | | | | | | | | | | | | | |
| | Apolipoprotein AI | | | | | | | | | | | | | | | | | |
| | Apolipoprotein B | | | | | | | | | | | | | | | | | |
| | aPTT | | | | | | | | | X | | | | | | | | |
| | AST (ASAT) | | | | | | | | | | | | | X | | | | |
| B | β-2-Microglobulin | | | | | | | | | | | | | | | | X | |
| | Benzoylcegonine | | | | | | | | | | | | | | | | | |
| | Bicarbonate | | | X | | | | | | | | | X | | | | | |
| | Bile Acids | | | | | | | | | | | | X | | | | | |
| | Bilirubin (Direct) | | | | | | | | | | | | X | | | | | |
| | Bilirubin (Total) | | | | | | | | | | | | X | | | | | |
| | Blood | | | | | | | | | | | | | | | | | |
| | BNP | | | | X | | | | | | | | | | | | | |
| | Buprenorphine | | | | | | | | | | | | | | | | | |
| C | CA15-3 | | | | | | | | | | | | | | | | X | |
| | CA19-9 | | | | | | | | | | | | | | | | X | |
| | CA125 | | | | | | | | | | | | | | | | X | |
| | Cadmium | | | | | | | | | | | | | | | | | |
| | Caffeine | | | | | | | | | | | | | | | | | |
| | Calcitonin | | | | | | | | | | | | | | | | | |
| | Calcium | | | | | | | | | | | | X | | X | | | |
| | Calcium, Adjusted* | | | | | | | | | | | | X | | | | | |
| | Calcium (Ionised) | | | X | | | | | | | | | X | | | | | |

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunosassay Speciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | |
|-----------------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|-----------------------------|---|
| | | | | | | | X | | | | | | | | | Anti-HSV-1 & 2 IgM Combined | A |
| | | | | | | | X | | | | | | | | | Anti-HSV1 IgG | |
| | | | | | | | X | | | | | | | | | Anti-HSV1 IgM | |
| | | | | | | | X | | | | | | | | | Anti-HSV2 IgG | |
| | | | | | | | X | | | | | | | | | Anti-HSV2 IgM | |
| | | | | X | | | | | | | | | | | | Anti-HTLV-I & 2 Combined | |
| | | | | X | | | | | | | | | | | | Anti-HTLV-I | |
| | | | | X | | | | | | | | | | | | Anti-HTLV-II | |
| | | | | | | | X | | | | | | | | | Anti-Measles IgG* | |
| | | | | | | | X | | | | | | | | | Anti-Mumps IgG* | |
| | | | | | | | X | | | | | | | | | Anti-Rubella IgG | |
| | | | | | | | X | | | | | | | | | Anti-Rubella IgM | |
| | | | | | | | | | | | | | | | | Anti-TG | |
| | | | | | | | X | | | | | | | | | Anti-VZV IgG* | |
| | | | | | | | | X | | | | | | | | Anti-thrombin III | |
| | | | | | | | X | | | | | | | | | Anti-Toxoplasma IgG | |
| | | | | | | | X | | | | | | | | | Anti-Toxoplasma IgM | |
| | | | | | | | | | | | | | | | | Anti-TPO | |
| | | | | | | | | | | | | | | | | Anti-TSH Receptor (TRAb) | |
| | X | | | | | | | | | | | | | | | Apolipoprotein AI | |
| | X | | | | | | | | | | | | | | | Apolipoprotein B | |
| | | | | | | | | | | | | | | | | aPTT | |
| | | | | | | | | | | | | | | | | AST (ASAT) | |
| | | | | | | | X | | | | | | | | | β-2-Microglobulin | B |
| | | | | | | | | | | | | | | X | | Benzoylcegonine | |
| | | | | | | | | | | | | | | | | Bicarbonate | |
| | | | | | | | | | | | | | | | | Bile Acids | |
| | | | | | | | | | | | | | | | | Bilirubin (Direct) | |
| | | | | | | | | | | | | | X | | | Bilirubin (Total) | |
| | | | | | | | | | | | | | X | | | Blood | |
| | | | | | | | | | | | | | | | | BNP | |
| | | | | | | | | | | | | | | X | | Buprenorphine | C |
| | | | | | | | | | | | | | | | | CA15-3 | |
| | | | | | | | | | | | | | | | | CA19-9 | |
| | | | | | | | | | | | | | | | | CA125 | |
| | | | | | | | | | | | | X | | | | Cadmium | |
| | | | | | | | | | X | | | | | | | Caffeine | |
| X | | | | | | | | | | | | | | | | Calcitonin | |
| | | | | | | | | | | | | | | | | Calcium | |
| | | | | | | | | | | | | | | | | Calcium, Adjusted* | |
| | | | | | | | | | | | | | | | | Calcium (Ionised) | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| | | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-1 + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality + |
|------------------------|--|---------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|
| C | Cannabinoids (THC) | | | | | | | | | | | | | | | | | |
| | Carbamazepine | | | | | | | | | | | | | | | | | X |
| | Carboxyhaemoglobin (COHb / HbCO) | | | | | | | | | X | | | | | | | | |
| | CEA | | | | | | | | | | | | | | | | | X |
| | Ceruloplasmin | | | | | | | | | | | | | | | | | |
| | Chloride | | | X | | | | X | | | | | X | | | X | | |
| | Cholesterol (Total) | | | | | | | | | | | | X | | | | | |
| | Cholinesterase | | | | | | | | | | | | X | | | | | |
| | Chromium | | | | | | | | | | | | | | | | | |
| | Ciclosporin | | | | | | | | | | | | | | | | | |
| | CK, Total | | | | | | X | X | | | | | X | | | | | |
| | CK-MB (Activity) | | | | | | X | X | | | | | | | | | | |
| | CK-MB (Mass) | | | | | | X | X | | | | | | | | | | |
| | CO ₂ , Total | | | X | | | | | | | | | | | | | | |
| | Cobalt | | | | | | | | | | | | | | | | | |
| | Complement C ₃ | | | | | | | | | | | | | | | | | |
| | Complement C ₄ | | | | | | | | | | | | | | | | | |
| | Conductivity | | | | | | | | | | | | | | | | | |
| | Copper | | | | | | | | | | | | | X | | | X | |
| | Cortisol | | | | | | | | | | | | | | | X | X | |
| | Cotinine | | | | | | | | | | | | | | | | | |
| | C-Peptide | | | | | | | | | | | | | | | | | X |
| | C-Reactive Protein (CRP) | | | | | | | | | | | | | | | | | X |
| | Creatinine | | | | | | | | | | | | | X | | X | | |
| | CYFRA 21-1 (Cytokeratin 19) | | | | | | | | | | | X | | | | | | |
| | D | D-3-Hydroxybutyrate | | | | | | | | | | | | X | | | | |
| | | d-Amphetamine | | | | | | | | | | | | | | | | |
| D-Dimer* ^Δ | | | | | | | X | X | | | | | | | | | | |
| Deoxyhaemoglobin (HHb) | | | | | | | | | | X | | | | | | | | |
| DHEA Unconjugated | | | | | | | | | | | | | | | | | X | |
| DHEA-Sulphate | | | | | | | | | | | | | | | | | X | |
| Digoxin | | | | | | | X | | | | | | | | | | X | |
| d-Methamphetamine | | | | | | | | | | | | | | | | | | |
| Dopamine | | | | | | | | | | | | | | | | X | | |
| E | | EDDP | | | | | | | | | | | | | | | | |
| | eGFR (estimated glomerular filtration rate)* | | | | | | | | | | | | X | | | | | |
| | Epinephrine | | | | | | | | | | | | | | | X | | |
| | ESR | | | | | | | | | | | X | | | | | | |
| | Ethanol | X | | | | | | | | | | | | | | | | |
| | Ethosuximide | | | | | | | | | | | | | | | | | |
| | Everolimus | | | | | | | | | | | | | | | | | |

^Δ Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunospeciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | |
|----------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|--|---|
| | | | | | | | | | | X | | | | | X | Cannabinoids (THC) | C |
| | | | | | | | | | | | | | | | | Carbamazepine | |
| | | | | | | | | | | | | | | | | Carboxyhaemoglobin (COHb / HbCO) | |
| | | | | | | | | | | | | | | | | CEA | |
| | | | | | | | X | | | | | | | | | Ceruloplasmin | |
| | | | | | | | | X | | | | | | | | Chloride | |
| | X | | | | | | | | | | | | | | | Cholesterol (Total) | |
| | | | | | | | | | | | | | | | | Cholinesterase | |
| | | | | | | | | | | | X | X | | | | Chromium | |
| X | | | | | | | | | | X | | | | | | Ciclosporin | |
| | | | | | | | | | | | | | | | | CK, Total | |
| | | | | | | | | | | | | | | | | CK-MB (Activity) | |
| | | | | | | | | | | | | | | | | CK-MB (Mass) | |
| | | | | | | | | | | | | | | | | CO ₂ , Total | |
| | | | | | | | | | | | X | X | | | | Cobalt | |
| | | | | | | | | X | | | | | | | | Complement C ₃ | |
| | | | | | | | | X | | | | | | | | Complement C ₄ | |
| | | | | | | | | | X | | | | | | | Conductivity | |
| | | | | | | | | | | X | X | X | | | | Copper | |
| | | | | | | | | | | | | | | | | Cortisol | |
| | | | | | | | | | | | | | | | X | Cotinine | |
| | | | | | | | | | | | | | | | | C-Peptide | |
| | | | | | | | X | | | | | | | | | C-Reactive Protein (CRP) | |
| | | | | | | | | | | | | | X | X | | Creatinine | |
| | | | | | | | | | | | | | | | | CYFRA 21-1 (Cytokeratin 19) | |
| | | | | | | | | | | | | | | | | D-3-Hydroxybutyrate | D |
| | | | | | | | | | | | | | | X | | d-Amphetamine | |
| | | | | | | | | | | | | | | | | D-Dimer* ^Δ | |
| | | | | | | | | | | | | | | | | Deoxyhaemoglobin (HHb) | |
| | | | | | | | | | | | | | | | | DHEA Unconjugated | |
| | | | | | | | | | | | | | | | | DHEA-Sulphate | |
| | | | | | | | | | | X | | | | | | Digoxin | |
| | | | | | | | | | | | | | | X | | d-Methamphetamine | |
| | | | | | | | | | | | | | | | | Dopamine | |
| | | | | | | | | | | | | | | X | | EDDP | |
| | | | | | | | | | | | | | | | | eGFR (estimated glomerular filtration rate)* | E |
| | | | | | | | | | | | | | | | | Epinephrine | |
| | | | | | | | | | | | | | | | | ESR | |
| | | | | | | | | | | | | | X | | | Ethanol | |
| | | | | | | | | | X | | | | | | | Ethosuximide | |
| X | | | | | | | | | | | | | | | | Everolimus | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^Δ Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| | | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-I + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality + |
|-----------|-------------------------|---------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|
| F | Factor II | | | | | | | | X | | | | | | | | | |
| | Factor IX | | | | | | | | X | | | | | | | | | |
| | Factor V | | | | | | | | X | | | | | | | | | |
| | Factor VII | | | | | | | | X | | | | | | | | | |
| | Factor VIII | | | | | | | | X | | | | | | | | | |
| | Factor X | | | | | | | | X | | | | | | | | | |
| | Factor XI | | | | | | | | X | | | | | | | | | |
| | Factor XII | | | | | | | | X | | | | | | | | | |
| | Ferritin | | | | | | | | | | | | | | | | | X |
| | Fibrinogen | | | | | | | | X | | | | | | | | | |
| | Folate | | | | | | | | | | | | | | | | | X |
| | Free Morphine | | | | | | | | | | | | | | | | | |
| | free β-hCG | | | | | | | | | | | | | | | | | |
| | Fructosamine | | | | | | | | | | | | X | | | | | |
| | FSH | | | | | | | | | | | | | | | | | X |
| G | γ-GT | | | | | | | | | | | | X | | | | | |
| | Galactose | | | | | | | | | | | | | | | | | |
| | Gastrin | | | | | | | | | | | | | | | | | |
| | Gentamicin | | | | | | | | | | | | | | | | | X |
| | Growth Hormone (GH) | | | | | | | | | | | | | | | | | X |
| | GLDH | | | | | | | | | | | | X | | | | | |
| | Glucose | | | X | | | X | | | | | | X | | | X | | |
| | | | | | | | | | | | | | | | | | | |
| H | Haematocrit (HCT) | | | | | | | | | | | | | | X | | | |
| | Haemoglobin (Hb) | | | | | | | | | | | | | | X | | | |
| | Total Haemoglobin (tHb) | | | | | | | | | X | | | X | | | | | |
| | Haptoglobin | | | | | | | | | | | | | | | | | |
| | HbA1c | | | | | | | | | | | | X | | | | | |
| | HBsAg | | | | | | | | | | | | | | | | | |
| | HBDH | | | | | | | | | | | | X | | | | | |
| | hCG | | | | | | | | | | | | | | | | | X |
| | HDL-Cholesterol | | | | | | | | | | | | | X | | | | |
| | Homocysteine | | | | | X | X | | | | | | | | | | | |
| | hsCRP | | | | | | X | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | |
| | I | IgA | | | | | | | | | | | | | | | | |
| IgE | | | | | | | | | | | | | | | | | | X |
| IGF-I | | | | | | | | | | | | | | | | | | X |
| IgG | | | | | | | | X | | | | | | | | | | |
| IgM | | | | | | | | | | | | | | | | | | |
| Inhibin A | | | | | | | | | | | | | | | | | | |
| Insulin | | | | | | | | | | | | | | | | | X | X |
| Iodine | | | | | | | | | | | | | | | | | | |

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunospeciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | | |
|----------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|--|-------------------------|---|
| | | | | | | | | | | | | | | | | | Factor II | F |
| | | | | | | | | | | | | | | | | | Factor IX | |
| | | | | | | | | | | | | | | | | | Factor V | |
| | | | | | | | | | | | | | | | | | Factor VII | |
| | | | | | | | | | | | | | | | | | Factor VIII | |
| | | | | | | | | | | | | | | | | | Factor X | |
| | | | | | | | | | | | | | | | | | Factor XI | |
| | | | | | | | | | | | | | | | | | Factor XII | |
| | | | | | | | | X | | | | | | | | | Ferritin | |
| | | | | | | | | | | | | | | | | | Fibrinogen | |
| | | | | | | | | | | | | | | | | | Folate | |
| | | | | | | | | | | | | | | | X | | Free Morphine | |
| | | | X | | | | | | | | | | | | | | free β-hCG | |
| | | | | | | | | | | | | | | | | | Fructosamine | |
| | | | | | | | | | | | | | | | | | FSH | |
| | | | | | | | | | | | | | | | | | γ-GT | G |
| X | | | | | | | | | | | | | | X | | | Galactose | |
| | | | | | | | | | | | | | | | | | Gastrin | |
| | | | | | | | | | | X | | | | | | | Gentamicin | |
| | | | | | | | | | | | | | | | | | Growth Hormone (GH) | |
| | | | | | | | | | | | | | | | | | GLDH | |
| | | | | | | | | | | | | | | | X | | Glucose | |
| | | | | | | | | | | | | | | | | | Haematocrit (HCT) | |
| | | | | | | | | | | | | | | | | | Haemoglobin (Hb) | |
| | | | | | | | | | | | | | | | | | Total Haemoglobin (tHb) | |
| | | | | | | | | X | | | | | | | | | Haptoglobin | |
| | | | | | | | | | | | | | | | | | HbA1c | |
| | | | | | X | | | | | | | | | | | | HBsAg | |
| | | | | | | | | | | | | | | | | | HBDH | |
| | | | | | | | | | | | | | | X | | | hCG | |
| | | X | | | | | | | | | | | | | | | HDL-Cholesterol | |
| | | | | | | | | | | | | | | | | | Homocysteine | |
| | | | | | | | | | | | | | | | | | hsCRP | |
| | | | | | | | | X | | | | | | | | | IgA | I |
| | | | | | | | | X | | | | | | | | | IgE | |
| | | | | | | | | | | | | | | | | | IGF- I | |
| | | | | | | | | X | | | | | | | | | IgG | |
| | | | | | | | | X | | | | | | | | | IgM | |
| | | | X | | | | | | | | | | | | | | Inhibin A | |
| | | | | | | | | | | | | | | | | | Insulin | |
| | | | | | | | | | | | X | X | X | | | | Iodine | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| | | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-I + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality + | |
|----------------------|--|---------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|--|
| I | Iron | | | | | | | | | | | | X | | | | | | |
| K | Kappa Light Chain (Free) | | | | | | | | | | | | | | | | | | |
| | Kappa Light Chain (Total) | | | | | | | | | | | | | | | | | | |
| | Ketones | | | | | | | | | | | | | | | | | | |
| L | Lactate | | | X | | | X | | | | | | X | | | | | | |
| | Lambda Light Chain (Free) | | | | | | | | | | | | | | | | | | |
| | Lambda Light Chain (Total) | | | | | | | | | | | | | | | | | | |
| | LD (LDH) | | | | | | | | | | | | X | | | | | | |
| | LDL-Cholesterol* ^Δ | | | | | | | | | | | | X | | | | | | |
| | Lead | | | | | | | | | | | | | | | | | | |
| | Leucocytes | | | | | | | | | | | | | | | | | | |
| | Lipase | | | | | | | | | | | | | X | | | | | |
| | Lipoprotein (a) | | | | | | | | | | | | | | | | | | |
| | Lithium | | | | | | | | | | | | | X | | | | | |
| | Lorazepam | | | | | | | | | | | | | | | | | | |
| | LSD | | | | | | | | | | | | | | | | | | |
| | Luteinising Hormone (LH) | | | | | | | | | | | | | | | | | X | |
| M | Magnesium | | | | | | | | | | | | X | | | X | | | |
| | Manganese | | | | | | | | | | | | | | | | | | |
| | MDMA | | | | | | | | | | | | | | | | | | |
| | Mean Cell Haemoglobin (MCH) | | | | | | | | | | | | | | X | | | | |
| | Mean Cell Haemoglobin Concentration (MCHC) | | | | | | | | | | | | | | X | | | | |
| | Mean Cell Volume (MCV) | | | | | | | | | | | | | | X | | | | |
| | Mean Platelet Volume (MPV) | | | | | | | | | | | | | | X | | | | |
| | Metanephrine | | | | | | | | | | | | | | | X | | | |
| | Methadone | | | | | | | | | | | | | | | | | | |
| | Methaemoglobin (MetHb) | | | | | | | | | X | | | | | | | | | |
| | Methotrexate | | | | | | | | | | | | | | | | | | |
| | Molybdenum | | | | | | | | | | | | | | | | | | |
| | Myoglobin | | | | | X | X | | | | | | | | | | | | |
| | N | NEFA | | | | | | | | | | | | X | | | | | |
| | | Nickel | | | | | | | | | | | | | | | | | |
| Nitrite | | | | | | | | | | | | | | | | | | | |
| Non-HDL Cholesterol* | | | | | | | | | | | | | X | | | | | | |
| Norepinephrine | | | | | | | | | | | | | | | | X | | | |
| Normetanephrine | | | | | | | | | | | | | | | | X | | | |
| Norpropoxyphene | | | | | | | | | | | | | | | | | | | |
| Nortriptyline | | | | | | | | | | | | | | | | | | | |
| NTproBNP | | | | | | | X | | | | | | | | | | | | |
| O | Oestradiol | | | | | | | | | | | | | | | | | X | |
| | Osmolality | | | | | | | | | | | | X | | X | | | | |

^Δ Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunoassay Speciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | |
|----------------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|--|---|
| | | | | | | | | | | | | | | | | Iron | I |
| | | | | | | | | X | | | | | | | | Kappa Light Chain (Free) | K |
| | | | | | | | | X | | | | | | | | Kappa Light Chain (Total) | K |
| | | | | | | | | | | | | | | X | | Ketones | |
| | | | | | | | | | | | | | | | | Lactate | L |
| | | | | | | | | X | | | | | | | | Lambda Light Chain (Free) | L |
| | | | | | | | | X | | | | | | | | Lambda Light Chain (Total) | L |
| | | | | | | | | | | | | | | | | LD (LDH) | |
| | | X | | | | | | | | | | | | | | LDL-Cholesterol* ^Δ | |
| | | | | | | | | | | | X | X | X | | | Lead | |
| | | | | | | | | | | | | | | X | | Leucocytes | |
| | | | | | | | | | | | | | | | | Lipase | |
| | | X | | | | | | | | | | | | | | Lipoprotein (a) | |
| | | | | | | | | | X | | | | | | | Lithium | |
| | | | | | | | | | | | | | | | X | Lorazepam | |
| | | | | | | | | | | | | | | | X | LSD | |
| | | | | | | | | | | | | | | | | Luteinising Hormone (LH) | |
| | | | | | | | | | | | X | | X | | | Magnesium | M |
| | | | | | | | | | | | X | X | X | | | Manganese | M |
| | | | | | | | | | | | | | | X | | MDMA | |
| | | | | | | | | | | | | | | | | Mean Cell Haemoglobin (MCH) | |
| | | | | | | | | | | | | | | | | Mean Cell Haemoglobin Concentration (MCHC) | |
| | | | | | | | | | | | | | | | | Mean Cell Volume (MCV) | |
| | | | | | | | | | | | | | | | | Mean Platelet Volume (MPV) | |
| | | | | | | | | | | | | | | | | Metanephrine | |
| | | | | | | | | | | | | | | | X | Methadone | |
| | | | | | | | | | | | | | | | | Methaemoglobin (MetHb) | |
| | | | | | | | | | X | | | | | | | Methotrexate | |
| | | | | | | | | | | | | | X | | | Molybdenum | |
| | | | | | | | | | | | | | | | | Myoglobin | |
| | | | | | | | | | | | | | | | | NEFA | N |
| | | | | | | | | | | | | X | X | | | Nickel | N |
| | | | | | | | | | | | | | | X | | Nitrite | |
| | | | | | | | | | | | | | | | | Non-HDL Cholesterol* | |
| | | | | | | | | | | | | | | | | Norepinephrine | |
| | | | | | | | | | | | | | | | | Normetanephrine | |
| | | | | | | | | | | | | | | | X | Norpropoxyphene | |
| | | | | | | | | | | | | | | | X | Nortriptyline | |
| | | | | | | | | | | | | | | | | NTproBNP | |
| | | | | | | | | | | | | | | | | Oestradiol | O |
| | | | | | | | | | | | | | | | | Osmolality | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^Δ Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| | | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-I + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality + | |
|--------------|----------------------------------|-----------------------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|---|
| O | Osteocalcin | | | | | | | | | | | | | | | | | X | |
| | Oxalate | | | | | | | | | | | | | | | X | | | |
| | Oxazepam | | | | | | | | | | | | | | | | | | |
| | Oxygen Content (O2CT) | | | | | | | | | X | | | | | | | | | |
| | Oxygen Saturation (sO2 / Vol O2) | | | | | | | | | X | | | | | | | | | |
| | Oxyhaemoglobin (O2Hb / HbO2) | | | | | | | | | X | | | | | | | | | |
| P | P24* | | | | | | | | | | | | | | | | | | |
| | PAPP-A | | | | | | | | | | | | | | | | | | |
| | Paracetamol (Acetaminophen) | | | | | | | | | | | | | | | | X | | |
| | pCO ₂ | | | X | | | | | | | | | | | | | | | |
| | pH | | | X | | | | | | | | | | | | | | | |
| | Phencyclidine | | | | | | | | | | | | | | | | | | |
| | Phenobarbital | | | | | | | | | | | | | | | | | X | |
| | Phenytoin | | | | | | | | | | | | | | | | | X | |
| | Phosphate (Inorganic) | | | | | | | | | | | | X | | | X | | | |
| | Plasma Renin Activity | | | | | | | | | | | | | | | | | | |
| | Plasminogen | | | | | | | | X | | | | | | | | | | |
| | Plateletcrit (PCT) | | | | | | | | | | | | | | X | | | | |
| | Platelets (PLT) | | | | | | | | | | | | | | X | | | | |
| | pO ₂ | | | | X | | | | | | | | | | | | | | |
| | Potassium | | | X | | | | | | | | | | X | | | X | | |
| | Prealbumin (Transthyretin) | | | | | | | | | | | | | | | | | | |
| | Primidone | | | | | | | | | | | | | | | | | | |
| | Procalcitonin | | | | | | | | | | | | | | | | | | X |
| | Progesterone | | | | | | | | | | | | | | | | | X | |
| | Prolactin | | | | | | | | | | | | | | | | | X | |
| | Protein (Total) | | | | | | | X | | | | | | X | | X | | | |
| | Protein C | | | | | | | | X | | | | | | | | | | |
| | Protein S | | | | | | | | X | | | | | | | | | | |
| | PSA (Free) | | | | | | | | | | | | | | | | | X | |
| | PSA (Total) | | | | | | | | | | | | | X | | | | X | |
| | PT (Including INR) | | | | | | | | X | | | | | | | | | | |
| | PTH | | | | | | | | | | | | | | | | | X | X |
| | R | Red Blood Cell Count (RBC) | | | | | | | | | | | | | | X | | | |
| | | Red Cell Distribution Width (RDW) | | | | | | | | | | | | | | X | | | |
| | | Renin (Direct Concentration) | | | | | | | | | | | | | | | | | |
| | | Retinol Binding Protein | | | | | | | | | | | | | | | | | |
| | | Rheumatoid Factor | | | | | | | | | | | | | | | | | |
| | S | Salicylic Acid | | | | | | | | | | | | | | | | X | |
| Secobarbital | | | | | | | | | | | | | | | | | | | |
| Selenium | | | | | | | | | | | | | | | | | | | |

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunoassay Speciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | |
|----------------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|-----------------------------------|---|
| | | | | | | | | | | | | | | | | Osteocalcin | O |
| | | | | | | | | | | | | | | | | Oxalate | |
| | | | | | | | | | | | | | | | X | Oxazepam | |
| | | | | | | | | | | | | | | | | Oxygen Content (O2CT) | |
| | | | | | | | | | | | | | | | | Oxygen Saturation (sO2 / Vol O2) | |
| | | | | | | | | | | | | | | | | Oxyhaemoglobin (O2Hb / HbO2) | |
| | | | | | X | | | | | | | | | | | P24* | P |
| | | | X | | | | | | | | | | | | | PAPP-A | |
| | | | | | | | | | | X | | | | | | Paracetamol (Acetaminophen) | |
| | | | | | | | | | | | | | | | | pCO ₂ | |
| | | | | | | | | | | | | | | X | | pH | |
| | | | | | | | | | | | | | | | X | Phencyclidine | |
| | | | | | | | | | | X | | | | | X | Phenobarbital | |
| | | | | | | | | | | X | | | | | | Phenytoin | |
| | | | | | | | | | | | | | | | | Phosphate (Inorganic) | |
| X | | | | | | | | | | | | | | | | Plasma Renin Activity | |
| | | | | | | | | | | | | | | | | Plasminogen | |
| | | | | | | | | | | | | | | | | Plateletcrit (PCT) | |
| | | | | | | | | | | | | | | | | Platelets (PLT) | |
| | | | | | | | | | | | | | | | | pO ₂ | |
| | | | | | | | | X | | | | | | | | Potassium | |
| | | | | | | | | | | X | | | | | | Prealbumin (Transthyretin) | |
| | | | | | | | | | | | X | | | | | Primidone | |
| X | | | | | | | | | | | | | | | | Procalcitonin | |
| | | | | | | | | | | | | | | | | Progesterone | |
| | | | | | | | | | | | | | | | | Prolactin | |
| | | | | | | | | | | | | | | X | | Protein (Total) | |
| | | | | | | | | | | | | | | | | Protein C | |
| | | | | | | | | | | | | | | | | Protein S | |
| | | | | | | | | | | | | | | | | PSA (Free) | |
| | | | | | | | | | | | | | | | | PSA (Total) | |
| | | | | | | | | | | | | | | | | PT (Including INR) | |
| | | | | | | | | | | | | | | | | PTH | |
| | | | | | | | | | | | | | | | | Red Blood Cell Count (RBC) | R |
| | | | | | | | | | | | | | | | | Red Cell Distribution Width (RDW) | |
| X | | | | | | | | | | | | | | | | Renin (Direct Concentration) | |
| | | | | | | | | X | | | | | | | | Retinol Binding Protein | |
| | | | | | | | | X | | | | | | | | Rheumatoid Factor | |
| | | | | | | | | | | X | | | | | | Salicylic Acid | S |
| | | | | | | | | | | | | | | X | | Secobarbital | |
| | | | | | | | | | | | X | | | | | Selenium | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

| | | Ammonia / Ethanol + | Anti-TSH Receptor + | Blood Gas | BNP + | Cardiac | Cardiac Plus • | Cerebrospinal Fluid + | Coagulation | CO-Oximetry + | CYFRA 21-I + | ESR + | General Clinical Chemistry | HbA1c | Haematology | Human Urine | Immunoassay | Immunoassay Speciality + | |
|-----------------------|------------------------------------|---------------------|---------------------|-----------|-------|---------|----------------|-----------------------|-------------|---------------|--------------|-------|----------------------------|-------|-------------|-------------|-------------|----------------------------|--|
| S | SHBG | | | | | | | | | | | | | | | | | X | |
| | Sirolimus | | | | | | | | | | | | | | | | | | |
| | Sodium | | | X | | | | X | | | | | X | | | X | | | |
| | Specific Gravity | | | | | | | | | | | | | | | | | | |
| | Syphilis | | | | | | | | | | | | | | | | | | |
| T | T ₃ (Free) | | | | | | | | | | | | X | | | | | X | |
| | T ₃ (Total) | | | | | | | | | | | | X | | | | | X | |
| | T ₄ (Free) | | | | | | | | | | | | X | | | | | X | |
| | T ₄ (Total) | | | | | | | | | | | | X | | | | | X | |
| | Tacrolimus | | | | | | | | | | | | | | | | | | |
| | Testosterone (Free)* | | | | | | | | | | | | | | | | | X | |
| | Testosterone (Total) | | | | | | | | | | | | | | | | | X | |
| | Thallium | | | | | | | | | | | | | | | | | | |
| | Theophylline | | | | | | | | | | | | | | | | | X | |
| | Thyroglobulin | | | | | | | | | | | | | | | | | X | |
| | TIBC | | | | | | | | | | | | X | | | | | | |
| | Tobramycin | | | | | | | | | | | | | | | | | | |
| | Total hCG | | | | | | | | | | | | | | | | | | |
| | Transferrin | | | | | | | | | | | | | X | | | | | |
| | Triglycerides | | | | | | | | | | | | | X | | | | | |
| | Troponin I | | | | | | X | X | | | | | | | | | | | |
| | Troponin T | | | | | | X | X | | | | | | | | | | | |
| | TSH | | | | | | | | | | | | X | | | | | X | |
| | TT | | | | | | | | X | | | | | | | | | | |
| | U | UIBC | | | | | | | | | | | | X | | | | | |
| Unconjugated Oestriol | | | | | | | | | | | | | | | | | | | |
| Urea | | | | | | | | | | | | | X | | | X | | | |
| Uric Acid | | | | | | | | | | | | | X | | | X | | | |
| Urobilinogen | | | | | | | | | | | | | | | | | | | |
| V | Valproic Acid | | | | | | | | | | | | | | | | | X | |
| | Vancomycin | | | | | | | | | | | | | | | | | X | |
| | Vitamin B12 | | | | | | | | | | | | | | | | | X | |
| | VMA | | | | | | | | | | | | | | | X | | | |
| W | Total White Blood Cell Count (WBC) | | | | | | | | | | | | | X | | | | | |
| Z | Zinc | | | | | | | | | | | X | | | | | | | |

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

PARAMETER INDEX

| Immunoassay Speciality 2 + | Immunosuppressant + | Lipid | Maternal Screening | Serology (EBV) + | Serology (HIV / Hepatitis) + | Serology (Syphilis) + | Serology (ToRCH) + | Specific Proteins | Sweat Testing + | Therapeutic Drug | Trace Elements in Blood + | Trace Elements in Serum + | Trace Elements in Urine + | Urinalysis + | Urine Toxicology + | | |
|----------------------------|---------------------|-------|--------------------|------------------|------------------------------|-----------------------|--------------------|-------------------|-----------------|------------------|---------------------------|---------------------------|---------------------------|--------------|--------------------|------------------------------------|---|
| | X | | | | | | | | | | | | | | | SHBG | S |
| | | | | | | | | | | | | | | | | Sirrolimus | |
| | | | | | | | | | | | | | | | | Sodium | |
| | | | | | | | | | | | | | | X | | Specific Gravity | T |
| | | | | | X | | | | | | | | | | | Syphilis | |
| | | | | | | | | | | | | | | | | T ₃ (Free) | T |
| | | | | | | | | | | | | | | | | T ₃ (Total) | |
| | | | | | | | | | | | | | | | | T ₄ (Free) | |
| | X | | | | | | | | | | | | | | | T ₄ (Total) | |
| | | | | | | | | | | | | | | | | Tacrolimus | U |
| | | | | | | | | | | | | | | | | Testosterone (Free)* | |
| | | | | | | | | | | | | | | | | Testosterone (Total) | |
| | | | | | | | | | | | | | X | | | Thallium | |
| | | | | | | | | | | X | | | | | | Theophylline | |
| | | | | | | | | | | | | | | | | Thyroglobulin | |
| | | | | | | | | | | | | | | | | TIBC | |
| | | | | | | | | | | X | | | | | | Tobramycin | |
| | | | X | | | | | | | | | | | | | Total hCG | |
| | | | | | | | X | | | | | | | | | Transferrin | |
| | X | | | | | | | | | | | | | | | Triglycerides | |
| | | | | | | | | | | | | | | | | Troponin I | |
| | | | | | | | | | | | | | | | | Troponin T | |
| | | | | | | | | | | | | | | | | TSH | |
| | | | | | | | | | | | | | | | | TT | |
| | | | | | | | | | | | | | | | | UIBC | U |
| | | | X | | | | | | | | | | | | | Unconjugated Oestriol | |
| | | | | | | | | | | | | | | | | Urea | |
| | | | | | | | | | | | | | | | | Uric Acid | V |
| | | | | | | | | | | | | | | X | | Urobilinogen | |
| | | | | | | | | | | X | | | | | | Valproic Acid | |
| | | | | | | | | | | X | | | | | | Vancomycin | |
| | | | | | | | | | | | | | | | | Vitamin B12 | W |
| | | | | | | | | | | | | | | | | VMA | |
| | | | | | | | | | | | | | | | | Total White Blood Cell Count (WBC) | Z |
| | | | | | | | | | | | X | X | | | | Zinc | |

+ = Not accredited

• = Accreditation status pending

* = Pilot study ongoing

PURPLE = The only parameters available on RQ9135/a

^a Pilot status only in certain programmes. Please check pages 30-34 for more information.

RELATED PRODUCTS

ACUSERA True Third Party Quality Controls

As a world leading manufacturer of multi-analyte true third party controls, thousands of laboratories rely on Randox to accurately assess test system performance and ultimately empower them with the confidence required to release patient test results. With more than 400 analytes available, the number of individual controls required to cover your test menu is significantly reduced while simultaneously reducing costs, time and storage space. A choice of formats is available, including liquid or lyophilised, which ensures flexibility and suitability for laboratories of all sizes and budgets. Many features of the Acusera range can help you to meet ISO 15189:2012 requirements:

- Designed to react to the test system in the same manner as a patient sample, helping to reduce inconvenient shifts in QC results when reagent batch is changed and ultimately providing a true indication of laboratory performance.
- The presence of analytes at key decision levels ensures accurate instrument performance and eliminates the need for additional low/high controls at extra expense.
- Manufactured independently from any instrument, the Acusera range delivers unbiased performance assessment with any instrument or method, while eliminating the need for multiple instrument specific controls.

Product Portfolio

Antioxidants | Blood Gas | Cardiac Markers | Routine Chemistry | Coagulation | Haematology | Diabetes | Immunoassay
Immunology | Infectious Diseases (Serology) | Lipids | POCT | Therapeutic Drugs | Toxicology | Urine Chemistry



Uniquely combining more than 100 analytes conveniently in a single control, laboratories can significantly reduce costs and consolidate without compromising on quality. As true third party controls, unbiased performance assessment with any instrument or method is guaranteed.


ACUSERA 24•7 Interlaboratory Data Management

Designed for use with the Acusera range of third party controls, the Acusera 24•7 software helps laboratories monitor and interpret their QC data. Access to an impressive range of features, including interactive charts, the automatic calculation of Measurement Uncertainty & Sigma Metrics and live peer group data generated from our extensive database of laboratory participants, ensures Acusera 24•7 is the most comprehensive package available.

- Advanced statistical analysis with automatic calculation of performance metrics including; Sigma, UM, TE & %Bias.
- Instantly discover how you compare to your peers with peer group statistics updated live in real-time reducing time and money spent troubleshooting.
- Interactive charts allowing you to add events and multiple data sets for quick and easy performance monitoring.
- Automated data import with bi-directional connection to LIMS (eliminating manual data entry).

Software Features

Dashboard | Result History | Interactive Levey-Jennings Charts | Interactive Histogram Charts
 Performance Summary Charts | Statistical Analysis Report | Statistical Metrics Report | Audit Trail Report
 Uncertainty of Measurement Report | Exception Report | Peer Group Statistics | Acusera Advisor



 'The laboratory shall have a procedure to prevent the release of patient results in the event of quality control failure. When the quality controls rules are violated and indicate that examination results are likely to contain significant errors the results shall be rejected... Quality Control data shall be reviewed at regular intervals to detect trends in examination performance.'

ISO 15189:2012

RANDOX - A GLOBAL DIAGNOSTIC SOLUTIONS PROVIDER

Randox has been supplying laboratories worldwide with revolutionary diagnostic solutions for over 35 years. Our experience and expertise allow us to create a leading product portfolio of high quality diagnostic tools which offer reliable and rapid diagnosis. We believe that by providing laboratories with the right tools, we can improve healthcare worldwide.



RX SERIES

Renowned for quality and reliability, the RX series combines robust hardware and intuitive software with the world leading RX series test menu comprising an extensive range of high quality reagents including routine chemistries, specific proteins, lipids, therapeutic drugs, drugs of abuse, antioxidants and diabetes testing. The RX series offers excellence in patient care delivering unrivalled precision and accuracy for results you can trust, guaranteeing real cost savings through consolidation of routine and specialised tests onto one single platform.



REAGENTS

Randox offers an extensive range of third-party diagnostic reagents which are internationally recognised as being of the highest quality; producing accurate and precise results. At Randox, we re-invest significantly in R&D to ensure we meet the ever-changing needs of the laboratory. Consequently, Randox offer a range of novel and superior performance assays, including; sdLDL-C, Lipoprotein (a), H-FABP, Adiponectin, Copper and Zinc. Applications are available detailing instrument-specific settings for the convenient use of Randox Reagents on numerous clinical chemistry analysers.



EVIDENCE SERIES

In 2002, Randox invented the world's first, Biochip Array Technology, offering highly specific tests, coupled to the highly sensitive chemiluminescent detection, providing quantitative results instantly changing the landscape of diagnostic testing forever. The Randox Evidence Series of multi-analyte immunoanalysers provide an unrivalled increase in patient information per sample offering diagnostic, prognostic and predictive solutions across a variety of disease areas with a highly advanced clinical and toxicology immunoassay test menu including cardiac, diabetes, drugs of abuse, metabolic and renal markers.

Contact us for more information on any of our products and services:

HEADQUARTERS

Randox Laboratories Ltd, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom

☎ +44 (0) 28 9445 4399 ✉ marketing@randox.com 🌐 riqas.com

INTERNATIONAL OFFICES



AUSTRALIA
Randox (Australia) Pty Ltd.
Tel: +61 (0) 2 9615 4640



BRAZIL
Randox Brasil Ltda.
Tel: +55 11 5181-2024



CHINA
Randox Laboratories Ltd.
Tel: +86 021 6288 6240



CZECH REPUBLIC
Randox Laboratories S.R.O.
Tel: +420 2 1115 1661



FRANCE
Laboratoires Randox
Tel: +33 (0) 130 18 96 80



GERMANY
Randox Laboratories GmbH
Tel: +49 (0) 215 1937 0611



HONG KONG
Randox Laboratories Hong Kong Limited
Tel: +852 3595 0515



ITALY
Randox Laboratories Ltd.
Tel: +39 06 9896 8954



INDIA
Randox Laboratories India Pvt Ltd.
Tel: +91 80 2802 5000



POLAND
Randox Laboratories Polska Sp. z o.o.
Tel: +48 22 862 1080



PORTUGAL
Irlindox Laboratorios Quimica Analitica Ltda
Tel: +351 22 589 8320



PUERTO RICO
Clinical Diagnostics of Puerto Rico, LLC
Tel: +1 787 701 7000



REPUBLIC OF IRELAND
Randox Teoranta
Tel: +353 7495 22600



SLOVAKIA
Randox S.R.O.
Tel: +421 2 6381 3324



SOUTH AFRICA
Randox Laboratories SA (Pty) Ltd.
Tel: +27 (0) 11 312 3590



SOUTH KOREA
Randox Korea
Tel: +82 (0) 31 478 3121



SPAIN
Laboratorios Randox S.L.
Tel: +34 93 475 09 64



SWITZERLAND
Randox Laboratories Ltd. (Switzerland)
Tel: +41 41 810 48 89



UAE
Randox Medical Equipments Trading LLC
Tel: +971 55 474 9075



USA
Randox Laboratories-US, Ltd.
Tel: +1 304 728 2890



VIETNAM
Randox Laboratories Ltd. Vietnam
Tel: +84 (0) 8 3911 0904

For technical support contact:

mail@riqas.com

RIQAS 



LT033 OCT20