

RANDOX

VERASTAT

Excellence at the Point of Care



Accurate & Reliable Results in Minutes

Innovation at the Point of Care



CONTENTS

- 01** BENEFITS
- 03** TESTING PROCESS
- 04** THE CASSETTE
- 05** TEST MENU
- 07** OUR TECHNOLOGY
- 08** SPECIFICATIONS
- 09** ORDERING INFORMATION
- 11** RANDOX - A GLOBAL DIAGNOSTIC SOLUTIONS PROVIDER
- 12** CONTACT US

BENEFITS

The Randox VeraSTAT is a simple, accurate, fast, portable point of care device based on cathodic electrochemiluminescence technology. Designed to offer users the next generation of rapid diagnosis providing results in as little as 6 minutes. The superiority of the VeraSTAT allows for users to overcome performance limitations of previous generation tests relating to sensitivity, accuracy, ease of use and cost efficiency. Combined with a versatile test menu, the Randox VeraSTAT is built to outshine and enhance detection in any setting.



ACCURATE RESULTS IN AS LITTLE AS 6 MINUTES

Eliminate delays in sending samples to the lab and facilitate immediate decision making at the point of care.



SIMPLE AND EASY TO USE

Intuitive user interface guides the operator through the entire testing process.



SMALL FOOTPRINT ANALYSER EASILY TRANSPORTED TO MULTIPLE TESTING LOCATIONS

Lightweight, portable and convenient the Randox VeraSTAT can be used in a variety of locations delivering diagnostic results when and where they are required.



ALL-IN-ONE CASSETTE

All necessary reagents are conveniently included in each single use, sealed cassette with no preparation required.



VERSATILE TEST MENU

Offering the differentiation of bacterial and viral acute respiratory tract infections, with a wide range of upcoming immunoassay, protein, inflammatory and diabetes disease markers



EXPORT RESULTS

The Radox VeraSTAT allows for results to be exported via bluetooth connectivity.



SMALL SAMPLE VOLUME


Requires just a few drops of whole blood or serum, ideal for use at the point of care.




BUILT IN QC MODULE

Dedicated QC module allowing QC history to be easily accessed.

THE TESTING PROCESS



THE CASSETTE



TEST MENU



INFLAMMATION AND INFECTION

C Reactive Protein (CRP)

C-reactive protein (CRP) is the dominant acute phase protein often used to guide treatment of a bacterial infection or inflammation associated with tissue injury, inflammatory disorders and associated diseases.

CRP measurements can also be used to assess the risk of cardiovascular disease, infection in neonates and for the early detection of renal allograft rejection.

CRP combined with MxA on the Randox VeraSTAT offer laboratories a powerful method in the fight against the antimicrobial resistance problem globally.

Sample Type: Whole Blood

Sample Volume: 5 µl

Assay Time: 6 mins

Working Range: 5-200 mg/L

Storage: 2-8 °C

Ordering Information: VS1003

Myxovirus Resistance Protein 1 (MxA)

Protein (MxA) is an informative general marker for the most common acute viral infections. As a key mediator of infection, MxA protein has the potential to greatly enhance the rapid distinction between viral and bacterial respiratory infections.

Sample Type: Whole Blood

Sample Volume: 10 ul

Assay Time: 11 mins

Working Range: 100-1000 µg/L

Storage: 2-8 °C

Ordering Information: VSI004

TEST MENU



DIABETES MARKERS

Glycated Haemoglobin (HbA1c) ***In Development**

For the quantitative determination of Hemoglobin A1c in whole blood samples.



CARDIAC MARKERS

Creatine Kinase MB (CKMB) ***In Development**


For the quantitative measurement of CK-MB in whole blood samples to assist diagnoses of an acute myocardial infarction.



OUR TECHNOLOGY

The Radox VeraSTAT uses patented Cathodic Electrochemiluminescence (C-ECL)

1. The target analyte in the patient sample reacts with the labelled antibody and capture antibody.
2. After the reaction, unbound or excess labelled antibody is washed away.
3. The labelled antibody complex is excited with electricity. Resulting electrochemiluminescence is measured.



SPECIFICATIONS



Dimensions	Height: 16cm, Depth: 25cm, Width: 18.5cm
Weight	3.5kg
Power	100-240 V AC (50/60 Hz)
Interface	Colour touch screen NFC reader included
Operational Environment	Temperature: +18 °C to +35 °C Maximum Relative Humidity: 80% Location: Dry, Clean, leveled surface away from direct sunlight, mechanical and magnetic interference
Memory	Stores up to 1000 results
Connectivity	Bluetooth connectivity allows results to be exported
Methodology	Cathodic Electrochemiluminescence (C-ECL)
Sample Type	Whole Blood, Serum
Results	Quantitative and qualitative
Sample Volume	Minimum sample volume 5ul
Operating Temperature	+18 °C to +35 °C

ORDERING INFORMATION



With almost 40 years' experience in the diagnostics industry, we fully understand the demands and needs of each customer. We recognise the importance of maintaining a consistent workflow of high quality results with minimum downtime. Our global network of highly trained engineers and customer support scientists are on hand to work with you in preserving the continuity of your operations while maximising the potential of your Randox instrument.

ORDERING INFORMATION

Randox VeraSTAT

VS1000

For more information contact your local representative, or email us at: marketing@randox.com



RANDOX - A GLOBAL DIAGNOSTIC SOLUTIONS PROVIDER

Randox has been supplying laboratories worldwide with revolutionary diagnostic solutions for almost 40 years. Our experience and expertise allow us to create a leading product portfolio of high quality diagnostic tools which offer reliable and rapid diagnosis. We believe that by providing laboratories with the right tools, we can improve healthcare worldwide.

REAGENTS



Randox offers an extensive range of open channel diagnostic reagents which are internationally recognised as being of the highest quality; producing accurate and precise results. At Randox, we re-invest significantly in R&D to ensure we meet the ever-changing needs of the laboratory. Consequently, Randox offer a range of novel and superior performance assays, including: sdLDL-C, Lipoprotein (a), Copper and Zinc. Applications are available detailing instrument-specific settings for the convenient use of Randox Reagents on numerous clinical chemistry analysers.

INTERNAL QUALITY CONTROL



Acusera third party quality controls are made using the highest quality material ensuring they react like a real patient sample. With more than 390 analytes available across the Acusera range we can uniquely reduce the number of controls required while reducing costs and time. Our product range includes clinical chemistry, immunoassay, urine, immunology and more. Qnostics molecular controls for infectious disease testing are designed to meet the demand of today's molecular diagnostics laboratory while effectively monitoring the entire testing process. Our whole pathogen molecular controls comprise hundreds of characterised viral, bacterial and fungal targets.

EXTERNAL QUALITY ASSESSMENT



RIQAS is the world's largest international EQA scheme with more than 55,000 participants worldwide. 39 comprehensive, yet flexible programmes cover a wide range of clinical diagnostic testing including chemistry, immunoassay, cardiac, urine, serology and more. Our programmes benefit from a wide range of concentrations, frequent reporting, rapid feedback and user-friendly reports. QCMD, our range of EQA programmes for molecular infectious disease testing, features a whole pathogen matrix making them a true test of patient sample analysis. With access to over 90 programmes including blood borne viruses, respiratory diseases, multi-pathogen infections and more, there is something for every laboratory.

RX SERIES



Renowned for quality and reliability, the RX series combines robust hardware and intuitive software with the world leading RX series test menu comprising an extensive range of high quality reagents including routine chemistries, specific proteins, lipids, therapeutic drugs, drugs of abuse, antioxidants and diabetes testing. The RX series offers excellence in patient care delivering unrivalled precision and accuracy for results you can trust, guaranteeing real cost savings through consolidation of routine and specialised tests onto one single platform.

CONTACT US

Contact us for more information on any of our products and services:

HEADQUARTERS

Randox Laboratories Ltd, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom
+44 (0) 28 9442 2413 • marketing@randox.com • randox.com

INTERNATIONAL OFFICES



AUSTRALIA
Randox (Australia) Pty Ltd.
Tel: +61 (0) 2 9615 4640



BRAZIL
Randox Brasil Ltda.
Tel: +55 11 5181 2024



CHINA
Randox Laboratories Ltd.
Tel: +86 021 6288 6240



CZECH REPUBLIC
Randox Laboratories S.R.O.
Tel: +420 2 1115 1661



FRANCE
Laboratoires Randox
Tel: +33 (0) 130 18 96 80



GERMANY
Randox Laboratories GmbH
Tel: +49 (0) 215 1937 0611



HONG KONG
Randox Laboratories Hong Kong Limited
Tel: +852 3595 0515



ITALY
Randox Laboratories Ltd.
Tel: +39 06 9896 8954



INDIA
Randox Laboratories India Pvt Ltd.
Tel: +91 80 67515000



POLAND
Randox Laboratories Polska Sp. z o.o.
Tel: +48 22 862 1080



PORTUGAL
Irlindox Laboratorios Quimica Analitica Ltda
Tel: +351 22 589 8320



PUERTO RICO
Clinical Diagnostics of Puerto Rico, LLC
Tel: +1 787 701 7000



REPUBLIC OF IRELAND
Randox Teoranta
Tel: +353 7495 22600



SLOVAKIA
Randox S.R.O.
Tel: +421 2 6381 3324



SOUTH AFRICA
Randox Laboratories SA (Pty) Ltd.
Tel: +27 (0) 11 312 3590



SOUTH KOREA
Randox Korea
Tel: +82 (0) 31 478 3121



SPAIN
Laboratorios Randox S.L.
Tel: +34 93 475 09 64



SWITZERLAND
Randox Laboratories Ltd. (Switzerland)
Tel: +41 41 810 48 89



UAE
Randox Medical Equipments Trading LLC
Tel: +971 55 474 9075



USA
Randox Laboratories-US, Ltd.
Tel: +1 304 728 2890



VIETNAM
Randox Laboratories Ltd. Vietnam
Tel: +84 (0) 8 3911 0904

FOR TECHNICAL SUPPORT CONTACT:

technical.services@randox.com

RANDOX



Copyright © Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151 6827 08. In EU VAT number: XI 151 6827 08. Product availability may vary from country to country. Some products may be for Research Use Only. For more information on product application and availability, please contact your local Randox Representative.

LT809 DEC23