

REGISTRATION INSTRUCTIONS & RIQAS POLICIES

CRITERIA FOR PARTICIPATION

This programme is available to any laboratory running the Serology (SYPHILIS) parameters listed in this document. Quantitative results will be accepted on this programme. Titre results will also be accepted. **RATIO results will not be accepted.**

INTRODUCTION

Method questionnaires are available for all routine RIQAS Programmes and are reviewed and updated every month, as indicated by the issue date at the bottom of every page. They are designed to allow you to register for this RIQAS Programme and to inform you of RIQAS protocols and policies. It is important that you read and understand all the information in these introductory pages before completing the enrolment document, which forms the basis of your registration and contract with RIQAS. If you have any questions or concerns about any of the information presented in this document, please contact RIQAS either directly or through your local Randox Laboratories representative. RIQAS Calendar dates and information about the RIQAS portfolio of products can be found at www.randox.com/external-quality-assessment.

REGISTRATION INSTRUCTIONS

NOTE: IF A REGISTERED PARTICIPANT DOES NOT PARTICIPATE FOR A CYCLE, THEY WILL BE EXPECTED TO COMPLETE NEW ENROLMENT DOCUMENTS IN ORDER TO RE-JOIN THE PROGRAMME.

METHOD QUESTIONNAIRE:- To be retained by participant

This method questionnaire should be completed and retained by you for your records. Please ensure that you complete the method questionnaire in full. Your details will help us to classify your results correctly and thus provide you with useful statistical data.

In order to fully complete this questionnaire you will also need a copy of the RIQAS Instruments and Reagent Suppliers which is available to download from the Randox website (www.randox.com/external-quality-assessment). Please ensure you have this list available when completing this questionnaire.

Following this introduction section is the method questionnaire which indicates the method codes available for each parameter along with the standard RIQAS unit. On the method questionnaire for each parameter you wish to run please tick the method appropriate to you, then state your instrument code, reagent code, and the units that you use in your laboratory if they are different from the RIQAS standard units. If codes are not available for your assay, please state the details of your method clearly in the section at the end of the enrolment document.

Once your method questionnaire has been completed, you must transfer the information onto your enrolment document.

ENROLMENT DOCUMENT:- To be returned to RIQAS

Please be aware that it may take up to 3 weeks to process enrolment documents if you are not entering your own assay details. When registering RIQAS enrolment documents, it is recommended that you state business contact details, rather than personal.

A. LABORATORY REFERENCE NUMBER

On receipt of an enrolment document, each participant is assigned a **laboratory reference number** which consists of a **participant number** which is unique to your laboratory and a **registration letter** which is assigned for each new registration we receive from you. If you are a current or previous participant, please state your **participant number** on the enrolment document. If you do not have a Laboratory Reference Number, this will be generated by RIQAS when you register for the first time. Please quote this number on all correspondence with RIQAS.

B. CYCLE/PRODUCT REQUIREMENTS

Please tick the cycles you wish to subscribe for. If there is more than one kit/product offered for the programme, please also tick the kit you wish to subscribe for.

C. PRIMARY CONTACT DETAILS

It is important to state the full address details of the Quality Assessment Officer or contact person who will receive all correspondence during the cycle. Please also state the company name of the Randox representative who is supplying you with the RIQAS product under 'Randox Office/Distributor'

Please inform RIQAS of any change to contact details as soon as possible.

D. RIQASNet

RIQASNet is a web-based online method for result entry / method changes and additions of parameters / viewing of released reports. To access RIQASNet go to www.riqas.net. Internet access and login details are required for RIQASNet and Adobe Reader is required for viewing reports. Your initial login information and password will be supplied by RIQAS. Once you have logged in for the first time you will be able to change your RIQASNet password. If you forget your password please follow the 'Forgotten Password' link. Your login information will be based on the 1st email address you supply on your enrolment document. A PDF copy of the report will be sent to this address and can also be sent to 2 other email addresses. These addresses should be stated on your enrolment document.

E. PDF REPORTS

Reports are sent as PDF files. These files can be sent to up to 3 email addresses. Adobe Reader is required to view the reports. The email addresses to which reports are sent can be reviewed and changed on RIQASNet.

F. CUSTOMER DECLARATION

The declaration indicates that by submitting your enrolment document to RIQAS, either directly or via your local Randox representative, you have read and understood the RIQAS policies stated in the most recent Method Questionnaire associated with this programme. You understand that the submission of your enrolment document to RIQAS marks the beginning of an on-going agreement, and you will be automatically enrolled in subsequent cycles of this programme until we receive written confirmation of your cancellation. This should be received 12 weeks prior to the month in which the cycle starts. You understand that you must inform RIQAS of any changes to your contact details, assay details or contract status. You authorise Randox Laboratories Ltd. to send communication related to the products and service provided to the e-mail or postal addresses stated on your submitted enrolment document. You understand that you are permitted to request disclosure of, change or erase personal details held by Randox Laboratories Ltd. at any time. Note: Method questionnaires are updated every month and the issue date is stated on every questionnaire and enrolment document.

G. REGISTRATION OF ASSAY DETAILS

Labs can register their assay details using RIQASNet or can complete the 'Registration of Assay Details' section of the enrolment document. Labs should tick the appropriate box under the 'Registration of Assay Details' section of the enrolment document. If a lab wishes RIQAS to register their assay details, they should complete the Registration of Assay Details section using the codes from this method questionnaire and the Instrument/Reagent Supplier Book.

Once a participant has registered they will receive an email containing their RIQASNet login information. Once you have successfully logged in to RIQASNet you will see your various laboratory reference numbers for each registered programme. If you have opted to add parameters/assay details using RIQASNet, please do so as soon as possible (see below).

If no code is available for your assay, please state the details of your method clearly in the section at the end of the enrolment document or follow the instructions on RIQASNet.

For Ortho-Clinical Diagnostics VITROS registrations, please state the 2 digit slide Generation number for each analyte.

If units other than the standard RIQAS units are used, please specify these in the boxes supplied.

ONCE COMPLETED, THE ENROLMENT DOCUMENT SHOULD BE SENT TO RIQAS FOR REGISTRATION.

H. UPDATING ASSAY DETAILS

It is possible to change your unit, method, instrument or reagent classification during a cycle.

Method Changes via RIQASNet: These can be made in the Assay Details section of the Data Entry menu. A list of your registered laboratory reference numbers will appear on screen. Select the laboratory reference number for which you would like to change the assay details. A current list of assay details will appear, click on the appropriate parameter. To change the details click the arrow box on the appropriate details and select a new one. Save the changes and submit them to RIQAS. Changes will not be instantaneously updated on RIQASNet but will be uploaded onto RIQASNet usually within 3 working days. It is possible to submit results and method changes together as method changes will be made before results are entered in to the RIQAS database.

I. ADDITION OF PARAMETERS / ASSAY DETAILS

Adding Parameters via RIQASNet: Parameters can be added using the Assay Details section of the Data Entry menu. A list of your registered laboratory reference numbers will appear on screen. Select the laboratory reference number for which you would like to add the assay details. At the top of the screen is 'Add Parameter'. Click on this and a list of parameters you are not registered for will appear. Select the parameter you wish to add and click the arrow box on the appropriate details and select your assay details. Save the changes and submit them to RIQAS. As above, additions will be available on RIQASNet usually within 3 working days.

ORDERING RIQAS PRODUCTS

Please ensure your purchase order for each cycle is placed with your local Randox representative 12 weeks prior to the month in which the cycle starts. This will ensure sufficient time to process and despatch your kit(s) to you. Participants from UK or Ireland may order products directly from RIQAS with an official order number. Orders received within 12 weeks of the start of the cycle will be processed with an additional administration fee. Current prices of RIQAS products are available from your local Randox Laboratories representative.

It may be possible to order RIQAS products during a cycle, subject to availability. Please contact your local Randox representative for more information.

SHIPPING AND RECEIPT OF RIQAS PRODUCTS

Provided that you have ordered sufficiently in advance, your RIQAS kit(s) will be shipped to you to arrive before the analysis date of the first sample in the kit. If you do not receive your kit(s) before this time, please contact your local Randox representative.

On RIQASNet please access your account and download the relevant Instructions For Use (IFU) document for the programme and cycle purchased. The IFU includes material characteristics, preparation, stability, storage and safety information. On receipt of your RIQAS kit, please check that:

- a) it is the product you ordered
- b) the correct number of samples are present as indicated on the IFU
- c) the samples have the appearance as indicated on the IFU and that none of them are damaged

Please notify your local Randox representative immediately if any of these are incorrect.

Please ensure that the product is immediately stored according to the recommendations on the package labelling.

ASSAY OF SAMPLES & RETURN OF RESULTS

Carefully read the instructions stated on the Instructions for Use (IFU) prior to preparation and assay of RIQAS samples. **These are available on RIQASNet only.** The RIQAS samples should be assayed at the recommended time specified on the IFU. Following appropriate preparation, samples should be treated as routine, unless otherwise stated on the IFU. Please assay the samples on or before the recommended date for analysis and forward your results to RIQAS by no later than **17:00 GMT on the FINAL DATE**, as indicated in the IFU. Results are submitted via RIQASNet, which can be accessed once you have received log in details via email. This will include a link to RIQASNet Instructions for Use.

In keeping with the objectives of EQA schemes, participants should be aware that collusion and falsification of results is considered to be unethical and constitutes scientific fraud. RIQAS policies must ensure that a laboratory is unaware of RIQAS means for comparison before submitting their own results. Where a result is not submitted by the final date, a report will be issued, but the missing results will be indicated as "No return" or "N" throughout the RIQAS reports. RIQAS permits the submission of late or corrected results only under the circumstances described below. Requests for the submission of late or corrected results must be submitted in writing and in English on RIQAS Form No. 9277-RQ (either by the participant or their local Randox Representative) and must be approved by RIQAS Management. The form is available on www.riqas.net.

Requests for the submission of late results must be accompanied by evidence that an error has been made, and that the error has not been caused by the participant.

Requests for the correction or removal of erroneous results must be accompanied by evidence that the error was non-analytical, as defined on form 9277-RQ. RIQAS is obliged to inform country-specific regulatory bodies of requests for correction of results (if they request such information for laboratory monitoring purposes).

New reports will be re-issued for late or corrected results only where there has been an error made by Randox Laboratories HQ, Randox representatives or distributors.

LATE RESULTS

In general, late results will not be accepted after the final date.

Late results will only be accepted where there has been an error made by Randox Laboratories HQ, Randox representatives or distributors.

CORRECTED RESULTS

Laboratories may correct results only if it can be determined that the error was non-analytical and where the request for submission is within 4 weeks of the original final date. A laboratory may correct a result under the following circumstances:

- Reconstituting a sample in an incorrect volume before analysis
- Assaying and/or submitting the results for the wrong sample
- Making a transcription error - submission of an analyser print-out indicating that the analysis date was before the final date is required.

DESPATCH OF REPORTS

PDF reports will be emailed within 72 hours of the FINAL DATE and for those registered for RIQASNet the PDF reports will be available on RIQASNet shortly after.

USE OF RIQAS REPORTS

Participants have permission to make copies of their RIQAS reports for internal use and for regulatory purposes only. RIQAS reports must not be duplicated for external use without permission from the RIQAS Scheme Co-ordinator. Under no circumstances should information on RIQAS reports be taken out of context or falsified in any way. Information regarding the format of RIQAS Reports and the monitoring of EQA performance can be found in the RIQAS Brochure at www.randox.com/external-quality-assessment

CONFIDENTIALITY

Participation in any RIQAS programme is considered to be strictly confidential. Any data transfer or correspondence with participants, either directly or via local Randox representative, will be deemed confidential. Participants should be aware that regulatory authorities have the right to request an assessment of a participant's performance. Where regulatory authorities are to be provided with a participant's results, participants will be notified.

GENERAL DATA PROTECTION REGULATION 2018 & UK DATA PROTECTION ACT 2018

Randox Laboratories Ltd. complies with GDPR and the UK Data Protection Act and holds the minimum information required to maintain the contract with RIQAS customers. Contact details are required in order to effectively provide you with the RIQAS products and services. Participants are not under any obligation to provide personal information to enter into a contract with RIQAS. We recommend that business contact details are provided. All data associated with the provision of RIQAS is collated, stored and processed confidentially and securely, to avoid unlawful processing, accidental loss or damage.

CERTIFICATES OF PARTICIPATION

Complimentary certificates of participation for each RIQAS programme are made available on RIQASNet to participants at the **end of the current cycle**, provided that **at least 50%** of results have been returned. Participants who enrol mid-cycle will be eligible for a Certificate for Participation if they have participated in at least 50% of samples available for the remainder of the cycle since enrolment. The certificate will specify the cycle, programme and the LABORATORY / HOSPITAL NAME which is detailed in the certificate section of RIQASNet. At the end of a cycle, a list of all eligible labs will be exported from RIQASNet and certificates will be created according to these details. Please ensure all certificate details are up to date in your RIQASNet account.

PERFORMANCE SURVEILLANCE OF UK LABS

RIQAS is obligated to identify and report persistent poor performing UK labs to the National Quality Assessment Advisory Panel. Poor performers are identified as those failing to meet performance criteria agreed with NQAAP. The performance criteria is specified in all performance surveillance correspondence with participants, and is also available on request. Participants are initially informed of poor performance by letter. Failure to improve performance will prompt details to be forwarded to NQAAP. All information sent to participants and NQAAP is strictly confidential. Please contact RIQAS if you require further information on Performance Surveillance.

PARTICIPANT FEEDBACK, COMPLAINTS & APPEALS

In order to ensure that RIQAS provides an appropriate and satisfying service, participants are invited to complete a feedback survey on RIQASNet. You may contact us at any time during the cycle, should you have any requests for additional programmes or parameters or comments regarding existing programmes.

RIQAS makes every effort to ensure that the samples provided are clinically challenging to as many laboratory systems as possible. For details, please contact RIQAS either directly or through your local Randox representative.

Should the need arise, participants may raise requests or enquiries through correspondence with the local Randox Laboratories representative or by contacting RIQAS directly. Participants may appeal against the evaluation of their performance by completing a PARTICIPANT APPEALS FORM, 10770-RQ. Participants may raise a complaint in relation to the product or service provided by completing the PARTICIPANT COMPLAINTS FORM, 10772-RQ. These forms are available on RIQASNet, or on request from RIQAS.

SUB-CONTRACTING

RIQAS sub-contracts aspects of the scheme. RIQAS accepts responsibility for the sub-contractors' work and protocols are in place to ensure that sub-contractors are deemed competent.

OUR COMPETENCE AS A PROFICIENCY TESTING PROVIDER

On request, RIQAS is willing to co-operate with participants seeking evidence of our competence as a proficiency testing provider or information on the design and implementation of RIQAS Programmes.

DEVIATION FROM EXISTING POLICIES/SERVICE

If there is any deviation from the existing policies or service, participants will be notified either directly or via their local Randox representative.

COMMUNICATION

As part of the service provided by Randox Laboratories Ltd., participants may be contacted by e-mail regarding updates and new products, in line with Randox Laboratories Ltd. privacy policy, as stated in www.randox.com.

THIS PROGRAMME HAS NOT BEEN ACCREDITED TO ISO/IEC 17043:2010

Please contact RIQAS at

Tel: +44 (0) 28 9445 4399

E-Mail mail@riqas.com

RIQAS Scheme Co-ordinator: Sarah Fleck

RANDOX LABORATORIES LTD., 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom

RQ9154 - SEROLOGY (SYPHILIS) PROGRAMME

METHOD QUESTIONNAIRE

Syphilis - FTA-Abs U/ml

CODE	METHOD
VSAPD	<input type="checkbox"/> Alphadia IFA Anti-TP IgG assay
VSFBSI	<input type="checkbox"/> Biocientifica Inmunofluor FTA-Abs
VSFBSI	<input type="checkbox"/> Biomerieux FTA-Abs
VSFDBS	<input type="checkbox"/> Dasit MarBlot Syphilis FTA-Abs
VSFEBR	<input type="checkbox"/> Euroimmun FTA-Abs
VSFHEG	<input type="checkbox"/> Hemagen FTA-Abs
VSFMDR	<input type="checkbox"/> Mastafuor FTA-Abs
VSFODR	<input type="checkbox"/> Osrodek Diagnostyczny FTA-Abs
VSFVIM	<input type="checkbox"/> Viro Immun FTA-Abs
VSFZEU	<input type="checkbox"/> Zeus Scientific FTA-ABS IFA

Other methods, please specify on enrolment document

INSTRUMENT CODE

REAGENT CODE

OTHER UNITS, SPECIFY

Syphilis - Immunoassay U/ml

CODE	METHOD	CODE	METHOD
VSIARC	<input type="checkbox"/> Abbott Architect/ Alinity	VSIIEA	<input type="checkbox"/> Iason ELIASON ELISA
VSIACL	<input type="checkbox"/> Acon Laboratories Foresight EIA	VSIIIE	<input type="checkbox"/> InTec ELISA
VSIADA	<input type="checkbox"/> Adaltis ELAgen Syphilis Ab Kit	VSIL21	<input type="checkbox"/> Lab 21 EIA
VSIABC	<input type="checkbox"/> Autobio CLIA	VSILFT	<input type="checkbox"/> Lifotronic eCLIA
VSIBSS	<input type="checkbox"/> Beckman Coulter / Sentinel Syphilis TP Latex	VSIMIK	<input type="checkbox"/> Mikrogen Recomline Treponema
VSIBKE	<input type="checkbox"/> Biokit Bioelisa Syphilis	VSIMIN	<input type="checkbox"/> Mindray Anti-TP
VSIBLN	<input type="checkbox"/> Bioland Nanosign Syphilis	VSIODE	<input type="checkbox"/> Omega Diagnostics ELISA
VSIBMX	<input type="checkbox"/> Biomaxima Syphilis Test	VSIORE	<input type="checkbox"/> Orgenics ELISA
VSIBIT	<input type="checkbox"/> Biomerieux Trepanostika	VSIPME	<input type="checkbox"/> Plasmatec ELISA
VSIBRB	<input type="checkbox"/> Biorad Bioplex Syphilis test	VSIRSS	<input type="checkbox"/> Radim Syphilis Screening
VSIBRS	<input type="checkbox"/> Bio-Rad Syphilis Total Anti-body test	VSIC6	<input type="checkbox"/> Roche Cobas e601/602
VSICPE	<input type="checkbox"/> Captia ELISA	VSIE8	<input type="checkbox"/> Roche Cobas e402/e801
VSICCL	<input type="checkbox"/> Chemclin Immunoassay	VSIEYS	<input type="checkbox"/> Roche Elecsys / Cobas e411
VSICTE	<input type="checkbox"/> CTK Recombilisa Syphilis Ab ELISA	VSIRME	<input type="checkbox"/> Roche Modular
VSIDAI	<input type="checkbox"/> Diagnostic Automation Inc	VSIRST	<input type="checkbox"/> Roche/Sekisui TPLA
VSIDBC	<input type="checkbox"/> Diagnostic BioProbes CLIA	VSIRS2	<input type="checkbox"/> Roche/Sekisui TPLA2
VSIDPE	<input type="checkbox"/> Diagnostic BioProbes ELISA	VSISDE	<input type="checkbox"/> SD Syphilis ELISA 3.0
VSIDSE	<input type="checkbox"/> Diagnostics Systems ELISA	VSISAM	<input type="checkbox"/> Siemens Atellica IM
VSIDMT	<input type="checkbox"/> DiaMedix Trep-Sure ELISA	VSICEN	<input type="checkbox"/> Siemens Centaur
VSDIA	<input type="checkbox"/> DiaProph Med	VSISEZ	<input type="checkbox"/> Siemens Enzygnost
VSILIA	<input type="checkbox"/> Diasorin Liaison	VSIDP2	<input type="checkbox"/> Siemens Immulite 2000
VSILIX	<input type="checkbox"/> Diasorin Liaison XL	VSIIISAT	<input type="checkbox"/> Shenzhen I Flash Anti TP
VSIDTS	<input type="checkbox"/> Diasorin Treponema Screen	VSISHA	<input type="checkbox"/> Sysmex HISCL Anti TP Assay Kit
VSIDME	<input type="checkbox"/> Diasorin Murex	VSIVBE	<input type="checkbox"/> Vector-Best RecombiBest ELISA
VSIDCH	<input type="checkbox"/> Diesse Chorus	VSIVE	<input type="checkbox"/> Vircell ELISA
VSIDEZ	<input type="checkbox"/> Diesse Enzywell	VITROS	<input type="checkbox"/> Vitros Immunodiagnostic Syphilis kit
VSIEUE	<input type="checkbox"/> Euroimmun IFA	VITROTPA	<input type="checkbox"/> Vitros Syphilis TPA
VSIFLG	<input type="checkbox"/> Fujirebio Lumipulse G Series	VITROTEST	<input type="checkbox"/> Vitrotest ELISA
VSIGBE	<input type="checkbox"/> General Biologicals ELISA	WANTAI	<input type="checkbox"/> Wantai ELISA
VSIHUM	<input type="checkbox"/> Human ELISA	WIENER	<input type="checkbox"/> Wiener Lab ELISA
		VSIZAE	<input type="checkbox"/> ZAO Ecolab AgCL-RMP

Other methods, please specify on enrolment document

INSTRUMENT CODE

REAGENT CODE

OTHER UNITS, SPECIFY

RQ9154 - SEROLOGY (SYPHILIS) PROGRAMME

METHOD QUESTIONNAIRE

Syphilis - RPR U/ml

CODE	METHOD
VSRASR	<input type="checkbox"/> Abbott Syfacard-R
VSRAG	<input type="checkbox"/> Agappe RPR
VSART	<input type="checkbox"/> Arkray RPR Test
VSRASI	<input type="checkbox"/> ASI RPR screening test
VSRAMR	<input type="checkbox"/> Atlas Medical RPR Card Test
VSRAXR	<input type="checkbox"/> Axis-Shield RPR test
VSRBDM	<input type="checkbox"/> BD Macro-Vue RPR card test
VSRBBS	<input type="checkbox"/> Bio-Bacter Serobacter
VSRBKR	<input type="checkbox"/> Biokit Reditest RPR
VSRBLO	<input type="checkbox"/> Biolabo RPR-Carbon
VSRBMX	<input type="checkbox"/> BioMaxima RPR test
VSRBIM	<input type="checkbox"/> Biomerieux Nosticon II RPR
VSRBRS	<input type="checkbox"/> Bio-Rad RPR test card
VSRBRX	<input type="checkbox"/> Biorex Syphilis RPR
VSRBSS	<input type="checkbox"/> Biosystems RPR kit
VSRBTR	<input type="checkbox"/> Biotec RPR test kit
VSRCEG	<input type="checkbox"/> Cenogenics RPR test
VSRCHL	<input type="checkbox"/> Chemelex RPR Carbon
VSRCHR	<input type="checkbox"/> ChrnoLab RPR
VSRCMR	<input type="checkbox"/> Cormay RPR
VSRCRD	<input type="checkbox"/> Crescent Diagnostics RPR
VSRDAR	<input type="checkbox"/> Diagnostic Automation RPR kit
VSRDIA	<input type="checkbox"/> Dialab Syphilis RPR
VSDRD	<input type="checkbox"/> Diesse RPR DAT
VSRDDS	<input type="checkbox"/> DDS Diagnostics
VSRREKO	<input type="checkbox"/> EKOLab Syphilis RPR
VSRILT	<input type="checkbox"/> ELITech RPR
VSRIMR	<input type="checkbox"/> Emapol RPR kit
VSRFAD	<input type="checkbox"/> FAR Diagnostics RPR Carbon
VSRFDR	<input type="checkbox"/> Fortress Diagnostics/Biorex RPR Carbon
VSRFUR	<input type="checkbox"/> Fumouze RPR card test
VSRGID	<input type="checkbox"/> Gesse Diagnostics RPR Carbon
VSRGMT	<input type="checkbox"/> GM Test RPR
VSRHUM	<input type="checkbox"/> Human Syphilis RPR test
VSRIMR	<input type="checkbox"/> InterMedical RPR/Carbon
VSRWI	<input type="checkbox"/> Inverness Medical Wampole Impact RPR
VSRIBR	<input type="checkbox"/> Lab 21 RPR test kit
VSRILCH	<input type="checkbox"/> Linear Chemicals RPR Carbon
VSRILLR	<input type="checkbox"/> Lorne Laboratories RPR
VSRILTA	<input type="checkbox"/> LTA Syphilis RPR
VSRIMBR	<input type="checkbox"/> Mascia Brunelli RPR
VSRIMLT	<input type="checkbox"/> MonlabTest RPR
VSRINMN	<input type="checkbox"/> Nal von Minden NADAL RPR Test
VSRINP	<input type="checkbox"/> Niarmedic Plus RPR
VSRIMD	<input type="checkbox"/> Omega Diagnostics Immutrep (Semi-Quant)
VSRIPST	<input type="checkbox"/> Plasmatec RPR test kit
VSRIPSC	<input type="checkbox"/> Pulse Scientific Inc. RPR test
VSRIQCA	<input type="checkbox"/> Quimica Clinica Aplicada RPR
VSRIRGR	<input type="checkbox"/> RAL Gernon RPR
VSRIRBX	<input type="checkbox"/> Ranbaxy Trepostat RPR
VSRIRAN	<input type="checkbox"/> Randox RPR test
VSRIRLS	<input type="checkbox"/> Rapid Labs RPR
VSRIRD	<input type="checkbox"/> Reckon Diagnostics RPR test
VSRIDGT	<input type="checkbox"/> RFCL Diagnova Trepostat
VSRIRLP	<input type="checkbox"/> Rodelg Laboratorios RPR
VSRIRRP	<input type="checkbox"/> Roche/Sekisui RPR
VSRISBR	<input type="checkbox"/> Servibio RPR
VSRISGM	<input type="checkbox"/> SGM Italia RPR
VSRIKHB	<input type="checkbox"/> Shanghai Kehua Bio-tech RPR kit
VSRISDR	<input type="checkbox"/> Span Diagnostics RPR test
VSRISRC	<input type="checkbox"/> Spinreact RPR Carbon
VSRISRQ	<input type="checkbox"/> Stanbio RPR Quicktest
VSRITDT	<input type="checkbox"/> Teco Diagnostics RPR test
VSRITLC	<input type="checkbox"/> Tody Laboratories RPR Carbon
VSRITUL	<input type="checkbox"/> Tulip Carbogen
VSRIVBR	<input type="checkbox"/> Vector-Best RPR
VSRIVDC	<input type="checkbox"/> Vedalab RPR Carbon
VSRIWIE	<input type="checkbox"/> Wiener RPR slide test

Other methods, please specify on enrolment document

INSTRUMENT CODE	<input type="text"/>
REAGENT CODE	<input type="text"/>
OTHER UNITS, SPECIFY	<input type="text"/>

RQ9154 - SEROLOGY (SYPHILIS) PROGRAMME

METHOD QUESTIONNAIRE

Syphilis - Rapid Tests U/ml

CODE	METHOD
VSPAAD	<input type="checkbox"/> Abbott Determine Syphilis TP Test
VSPABO	<input type="checkbox"/> Abon Syphilis Ultra Rapid Test Device
VSPACS	<input type="checkbox"/> Acon Syphilis One Step
VSPACC	<input type="checkbox"/> Accurate One Step Syphilis
VSPACT	<input type="checkbox"/> Accurex Accutest Syphilis
VSPADB	<input type="checkbox"/> Ad-Bio Combo Rapid Test
VSPDET	<input type="checkbox"/> Alere/Inv Medical Determine Syphilis TP
VSPADS	<input type="checkbox"/> All Diag Syphilistop Optima
VSPABQ	<input type="checkbox"/> Amvi Biotech Syphilis Quick Test
VSPART	<input type="checkbox"/> Artron Syphilis Rapid Test
VSATM	<input type="checkbox"/> Atlas Medical Syphilis Ab Rapid Test
VSPAVC	<input type="checkbox"/> Avecon Syphilis Rapid test
VSPWAR	<input type="checkbox"/> Beijing Wantai Syphilis Rapid test
VSPBBT	<input type="checkbox"/> Biofocus Biotracer Syphilis Rapid Card
VSPBLN	<input type="checkbox"/> Bioland Nanosign Syphilis test strip
VSPBSU	<input type="checkbox"/> Bioplastic Syphilis Ultra Rapid Device
VSPBIP	<input type="checkbox"/> Bio-Rad ID-PaGIA
VSPBTR	<input type="checkbox"/> Biotest RightSign Syphilis Rapid Test
VSPBRT	<input type="checkbox"/> Biozek Syphilis Rapid Test Cassette
VSPCSR	<input type="checkbox"/> Citest Syphilis Rapid Test
VSPCTK	<input type="checkbox"/> CTK Biotech Rapid test
VSPDSU	<input type="checkbox"/> Diagnosis Syphilis Ultra Rapid
VSPDIP	<input type="checkbox"/> DiaMed ID PaGIA Test
VSPFDS	<input type="checkbox"/> Fortress Diagnostics Syphilis Strips
VSPHGS	<input type="checkbox"/> Healgen Scientific Syphilis Rapid Test
VSPHHS	<input type="checkbox"/> Human Hexagon Syphilis
VSPHRT	<input type="checkbox"/> Human Hexagon Syphilis Rapid test
VSPHU	<input type="checkbox"/> Humasis Syphilis Card
VSPHSU	<input type="checkbox"/> Hydrex Syphilis Ultra
VSPIBS	<input type="checkbox"/> IBS Syphilis Rapid test
VSPING	<input type="checkbox"/> Innogenetics INNO-LIA test strips
VSPINO	<input type="checkbox"/> Innovacon Instalert Rapid tests
VSPITO	<input type="checkbox"/> InTec One Step Syphilis test
VSPKLQ	<input type="checkbox"/> Koroglu Laboquick Anti Syphilis
VSPL21	<input type="checkbox"/> Lab 21 Rapid tests
VSPLRK	<input type="checkbox"/> Laborex Kenya Rapid Test
VSPMPA	<input type="checkbox"/> MP Diagnostics Assure Syphilis RT
VSPMDR	<input type="checkbox"/> Mikrogen Diagnostik Recline
VSPNMN	<input type="checkbox"/> Nal von Minden NADAL Syphilis test
VSPODI	<input type="checkbox"/> Omega Diagnostics Immunotrep (Slide)
VSPOMV	<input type="checkbox"/> Omega Diagnostics Visitect
VSPORR	<input type="checkbox"/> Orgenics Rapid tests
VSPQPD	<input type="checkbox"/> Qualpro Diagnostics Rapid test
VSPRLS	<input type="checkbox"/> Rapid Labs Syphilis
VSPRDS	<input type="checkbox"/> Reckon Diagnostics Strip test
VSPSDB	<input type="checkbox"/> SD Bioline Syphilis Rapid Test
VSPSDT	<input type="checkbox"/> Span Diagnostics Trust test
VSPTDO	<input type="checkbox"/> Teco Diagnostics One Step Cassette
VSPTYD	<input type="checkbox"/> Toyo Diagnostics Cassette
VSTRU	<input type="checkbox"/> Trueline Syphilis Ab Rapid Test
VSPVET	<input type="checkbox"/> VeroTest Syphilis Ab Rapid Test
VSPZSW	<input type="checkbox"/> Zephyr Syphicheck - WB

Other methods, please specify on enrolment document

INSTRUMENT CODE

REAGENT CODE

OTHER UNITS, SPECIFY

RQ9154 - SEROLOGY (SYPHILIS) PROGRAMME

METHOD QUESTIONNAIRE

Syphilis - TPHA U/ml

CODE	METHOD	CODE	METHOD
VSTAMT	<input type="checkbox"/> Atlas Medical TPHA kit	VSTSOB	<input type="checkbox"/> Sobioda TPHA
VSTAXD	<input type="checkbox"/> Axiom Diagnostic TPHA	VSTSRT	<input type="checkbox"/> Spinreact TPHA
VSTAXT	<input type="checkbox"/> Axis-Shield Microsoph TPHA	VSTTLT	<input type="checkbox"/> Tody Laboratories TPHA
VSTBIK	<input type="checkbox"/> Biokit Syphagen TPLH	VSTZAE	<input type="checkbox"/> ZAO Ecolab TPHA
VSTBLO	<input type="checkbox"/> Biolabo TPHA	VSTVLT	<input type="checkbox"/> Veda Lab TPHA Syphilis
VSTBIM	<input type="checkbox"/> Biomerieux TPHA		
VSTBRD	<input type="checkbox"/> Biorex Diagnostics TPHA		
VSTBRS	<input type="checkbox"/> Bio-Rad TPHA kit		
VSTBSS	<input type="checkbox"/> Biosystems TPHA kit		
VSTBTC	<input type="checkbox"/> Biotec TPHA kit		
VSTBPC	<input type="checkbox"/> BPC Biosed Syphilide TPHA		
VSTCYD	<input type="checkbox"/> Cypress Diagnostics TPHA		
VSTCHR	<input type="checkbox"/> Chorolab TPHA		
VSTDMS	<input type="checkbox"/> DAC SpectroMed TPHA-DAC		
VSTDIA	<input type="checkbox"/> Dialab Syphilis TPHA		
VSTDSF	<input type="checkbox"/> Diesse Syphilis Fast		
VSDTB	<input type="checkbox"/> Diesse TPHA DAT		
VSTEUR	<input type="checkbox"/> Eurobio TPHA test		
VSTELT	<input type="checkbox"/> ELITech TPHA		
VSTFAR	<input type="checkbox"/> FAR TPHA kit		
VSTFDT	<input type="checkbox"/> Fortress Diagnostics TPHA		
VSTFSP	<input type="checkbox"/> Fujirebio Serodia Particle Agg.		
VSTFST	<input type="checkbox"/> Futura Systems TPHA		
VSTHUM	<input type="checkbox"/> Human Syphilis TPHA liquid		
VSTIMT	<input type="checkbox"/> InterMedical Syphilis/TPHA		
VSTL21	<input type="checkbox"/> Lab 21 TPHA		
VSTLBK	<input type="checkbox"/> Labkit TPHA		
VSTLTA	<input type="checkbox"/> LTA Syphilis TPHA		
VSTLCH	<input type="checkbox"/> Linear Chemicals TPHA		
VSTLL	<input type="checkbox"/> Lorne Laboratories TPHA kit		
VSTMBT	<input type="checkbox"/> Mascia Brunelli TPHA		
VSTNMN	<input type="checkbox"/> Nal von Minden NADAL TPHA Test		
VSTNBM	<input type="checkbox"/> Newbio TPHA kit		
VSTOMD	<input type="checkbox"/> Omega Diagnostics Immutrep		
VSTOXO	<input type="checkbox"/> Oxoid TPHA test		
VSTPST	<input type="checkbox"/> Plasmatec TPHA test kit		
VSTRAN	<input type="checkbox"/> Randox TPHA kit		
VSTRLS	<input type="checkbox"/> Rapid Labs TPHA		
VSTSGM	<input type="checkbox"/> SGM Italia TPHA		
VSTSCG	<input type="checkbox"/> Siemens Cellognost		

Other methods, please specify on enrolment document

INSTRUMENT CODE

REAGENT CODE

OTHER UNITS, SPECIFY

RQ9154 - SEROLOGY (SYPHILIS) PROGRAMME

METHOD QUESTIONNAIRE

Syphilis - VDRL U/ml

CODE	METHOD
VSVADV	<input type="checkbox"/> All.Diaq VDRL check
VSVBCD	<input type="checkbox"/> Becton Dickinson VDRL
VSVBIO	<input type="checkbox"/> Bioclin Quibasa VDRL
VSVMD	<input type="checkbox"/> Biomed VDRL
VSVBIM	<input type="checkbox"/> Biomerieux VDRL
VSVBRS	<input type="checkbox"/> Bio-Rad VDRL Latex Pasteur
VSVBTC	<input type="checkbox"/> Biotec VDRL Kit
VSVBTV	<input type="checkbox"/> Britania VDRL - Sepssi
VSVCEG	<input type="checkbox"/> Cenogenics VDRL test
VSVCHR	<input type="checkbox"/> Chronolab VDRL
VSVDSY	<input type="checkbox"/> Diaqast Sypal CB
VSVDIA	<input type="checkbox"/> Dialab Syphilis VDRL
VSVDMA	<input type="checkbox"/> Diasorin Murex VDRL
VSVICV	<input type="checkbox"/> INCDMI Cantacuzino Spl VDRL
VSVIDV	<input type="checkbox"/> Investigacion Diagnosticos VDRL Test
VSVLBV	<input type="checkbox"/> Labtest VDRL
VSVLIC	<input type="checkbox"/> Licon VDRL-USR
VSVLL	<input type="checkbox"/> Lorne Laboratories Syphilis Serology Kit
VSVOMD	<input type="checkbox"/> Omega Diagnosticos Immutrep
VSVOXO	<input type="checkbox"/> Oxoid VDRL test
VSVPST	<input type="checkbox"/> Plasmatec VDRL test kit
VSVROD	<input type="checkbox"/> Rodelq Laboratories VDRL
VSVSDM	<input type="checkbox"/> Sclavo Diagnostics Microgen VDRL
VSVSCA	<input type="checkbox"/> Siemens VDRL Cardiolipin
VSVSDT	<input type="checkbox"/> Span Diagnostics TRUST VDRL
VSVSRV	<input type="checkbox"/> Spinreact VDRL
VSVTUG	<input type="checkbox"/> Tulip Carbogen
VSVTUT	<input type="checkbox"/> Tulip Trepolipin
VSVWIE	<input type="checkbox"/> Wiener Lab VDRL test

Other methods, please specify on enrolment document

INSTRUMENT CODE

REAGENT CODE

OTHER UNITS, SPECIFY
