


RANDOX

**RIQAS MICROBIOLOGY
(BACTERIAL IDENTIFICATION)
PROGRAMME**



RANDOX INTERNATIONAL QUALITY ASSESSMENT SCHEME

RIQAS MICROBIOLOGY (BACTERIAL IDENTIFICATION) PROGRAMME



Introducing our new RIQAS Microbiology programme. 3 samples will be distributed, twice per year with one sample being analysed every 2 months. Each strain is presented in a lyophilised pellet form within a self-contained tube containing a pre-filled ampoule for reconstitution, therefore helping to avoid cross contamination.

Case studies are provided with each sample that have been manufactured with robustly characterised strains.

Identification of the micro-organisms can be made at Gram positive/negative, Genus and Species level.

Multiple registrations mean that this programme can be used for training purposes to review different members of staff and their ability to correctly identify the organisms.

Antibiotic resistance assessment (Resistant/Intermediate/Susceptible) is provided for your own and other guidelines

Description	Frequency	Cat. No
Microbiology Programme	Bi-monthly, 3 samples x 2 distributions per year	RQ9197



Multiple registrations allow for comparison between lab staff



Submit results and view reports online via RIQAS.Net




Possible to be used as a training tool for new and current staff



Samples are manufactured to the highest quality.



Lyophilised material for enhanced stability



ACUSERA 24•7 Online QC software with real-time peer group statistics

Designed for use with the Acusera range of third party controls, the Acusera 24•7 software will help you monitor and interpret your QC data. Access to an impressive range of features, including interactive charts, the automatic calculation of Measurement Uncertainty & Sigma Metrics and live peer group data generated from our extensive database of laboratory participants, ensures Acusera 24•7 is the most comprehensive package available.



ACUSERA True third party controls offering complete test menu consolidation

Uniquely combining more than 100 analytes conveniently in a single control, laboratories can significantly reduce costs and consolidate without compromising on quality. As true third party controls, unbiased performance assessment with any instrument or method is guaranteed.

Randox Laboratories Ltd, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom
☎ +44 (0) 28 9442 2413 ✉ mail@riqas.com 🌐 randox.com



Copyright © 2022 Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151 6827 08. In EU VAT number: XI 151 6827 08. Product availability may vary from country to country. Some products may be for Research Use Only. For more information on product application and availability, please contact your local Randox Representative.